

Product Catalog

For Restorative Technologies



Important Product Information

Biomet 3i Restorative Products

For detailed product information on all Biomet 3i Restorative Products, refer to the Biomet 3i Restorative Products IFU (P-IIS086GR) located at ifu.biomet3i.com.

Table Of Contents

Introduction

Company Overview	i
Implant & Abutment System Connections	ii
BellaTek® Patient Specific Technologies Overview	iv
BellaTek Bars & Frameworks Design Options	v
Restorative Products Flowchart	vii

Certain Internal Connection:



External Hex Connection:



Certain Internal and External Hex Connections:



ICON KEY: These icons represent the connection types of the BIOMET **3i** Implant and Abutment Systems for both internal and external connection types represented throughout this catalog.

Restorative Products

Certain® Internal Connection:

Provisional Restorations

PreFormance® Posts	1
PreFormance Temporary Cylinders	2
Titanium Temporary Cylinders	2

Cement-Retained Restorations

GingiHue® Posts	3
GingiHue Post Components	4
UCLA Abutments	5
UCLA Abutment Components	6

Screw-Retained Restorations

UCLA Abutments	7
UCLA Abutment Components	8
Low Profile Abutments	9
Low Profile Abutment Components	10

Overdenture Restorations

LOCATOR® Abutments	11
LOCATOR Abutment Components	12

External Hex Connection:

Provisional Restorations

PreFormance Posts	13
PreFormance Temporary Cylinders	14
Titanium Temporary Cylinders	14

Cement-Retained Restorations

GingiHue Posts	15
GingiHue Post Components	16
UCLA Abutments	17
UCLA Abutment Components	18

Screw-Retained Restorations

UCLA Abutments	19
UCLA Abutment Components	20
Low Profile Abutments	21
Low Profile Abutment Components	22

Overdenture Restorations

LOCATOR Abutments	23
LOCATOR Abutment Components	24

Continued on next page. ➤

Table Of Contents (Cont'd)

Screws

Certain® Internal Connection:

Abutment, Retaining, Waxing & Try-In Screws	25
---	----

External Hex Connection:

Abutment, Retaining, Waxing & Try-In Screws	26
---	----

Miscellaneous Restorative Products

Mechanical Driver Kits	27
Drivers & Driver Tips	28
Navigator® System Instrumentation Overview	29
Tapered Navigator Laboratory Kit	30
Parallel Walled Navigator Laboratory Kit	30
Laboratory Tools	31
Miscellaneous Items	31

Patient Education Products 32

Additional Information

Index	33
Ordering Form	40
Ordering Information	42

Introduction

Company Overview



Providing Solutions – One Patient At A Time™

Mission Statement:

BIOMET **3i** will be recognized as a global leader in the oral reconstruction marketplace. We will achieve this through superior customer support and continuous product innovation.

Headquartered in Palm Beach Gardens, Florida, with operations throughout the world, BIOMET **3i** is one of the leading companies in the oral reconstruction market. Of equal importance, BIOMET **3i** brings the same innovation, high standards and comprehensive approach to customer service, dental practice and laboratory support and education.

BIOMET **3i** offers one of implant dentistry's most comprehensive lines of implants and abutments, augmented by a growing line of site preparation, digital, restorative and regenerative products.


BIOMET **3i** pioneered the development of biologically driven implants, winning worldwide acclaim for the microtextured OSSEOTITE® surface and then the nano-scaled, Bone Bonding® NanoTite™ Surface Implants. The tradition of implant innovation continues with the introduction of the **3i** T3® Implant, a contemporary hybrid that is designed for sustainable aesthetics.

This spirit of innovation continues with the introduction of DIEM®2, a protocol to enable clinicians to pursue immediate provisionalization options for full arch patient cases.



Introduction


Implant & Abutment System Connections




Certain® Internal Connection

The internal connection aspect of the implant and abutment are designed for ease of use and simplicity.

- QuickSeat® Connection: It Clicks! The Certain Implant and Abutment Systems feature the QuickSeat Connection. This unique connection produces an audible and tactile “click” that confirms placement of Impression Copings and Abutments.







- Abutment fingers cause the “click” and also provide retention for the prosthetic components in the implant before the screw is placed. A screw is needed to fully seat the components when the restoration is being tried in or definitively placed.



- The 6/12 hex inside the internal connection incorporates both a 6-point single and a 12-point double hex. The 6-point single hex has two functions: engaging the driver tip for mountless delivery during implant placement and providing anti-rotation for all straight abutments. The 12-point double hex provides 30° rotational positioning for pre-angled abutments.




- Impression Copings: Pick-Up and Twist Lock™ Transfer Copings are available in the Certain Internal Connection design, in three EP® (Emergence Profile) diameters and a straight diameter to match the corresponding Healing Abutment. The Impression Copings also “click” when properly seated. These copings have a different finger design that engages 2 mm internally into the implant, as compared to the abutment’s design that engages 4 mm deep into the implant. This shorter engagement into the implant is important because it allows off-axis draw for implant impressions that are divergent.

(purple)	(blue)	(yellow)	(green)
			
3.4 mm	4.1 mm	5 mm	6 mm

- Color-Coding: Certain Implants and all interfacing components are color-coded by restorative platform diameter for easy identification and selection.

Introduction


Implant & Abutment System Connections




■ External Hex Connection

The external hex connection implant and abutment designs have been widely used in implant dentistry.


- Gold Standard ZR (Zero Rotation) is a patented design incorporating machined microstops within the corners of the hex in the abutment. The design aims to reduce horizontal rotation between the implant and abutment. Gold Standard ZR is available on UCLA, GingiHue® and Conical Abutments.



- The 12-point double hex, which is machined in the 15° Pre-Angled GingiHue Abutment, provides rotational positioning in 30° increments on the implant hex.



- Precise Abutment Placement: The ASYST Abutment Placement Tool provides fingertip control for fast and easy abutment delivery. The patented packaging design makes abutment placement easier because the abutment is delivered sterile from the package directly to the implant site. The abutment and abutment screw are packaged inside the plastic ASYST seating device. The plastic seating device facilitates precise placement of the abutment on the implant, thus reducing chair time. The ASYST Tool is packaged with Low Profile, Standard and Conical Abutments.



- Impression Copings: Pick-Up and Twist Lock™ Transfer Copings are both available in the external connection design in three EP® (Emergence Profile) diameters and a straight diameter to match the corresponding healing abutment. The patented Twist Lock design provides a more accurate transfer of the implant connection as compared to other transfer (closed tray) techniques. This unique Twist Lock feature provides a series of undercuts to guide the coping into the impression, giving stability to the impression coping/analog prior to pouring the impression. As a result, the copings lock into orientation grooves upon clockwise rotation, thereby allowing for a tactile sensation of resistance that ensures an accurate hex transfer.

Introduction



Patient Specific Technologies Overview



BellaTek® Encode® Impression System

The BellaTek Encode Impression System eliminates the need for an implant-level impression while delivering a patient-specific restoration that should have the appropriate margin height and emergence contours for the patient. This accuracy is achieved from an impression of the BellaTek Encode Healing Abutment. Codes embedded on the occlusal surface of the BellaTek Encode Healing Abutment communicate the collar height, implant hex-orientation, platform diameter and interface (the Certain® Internal Connection or External Hex Connection) of the implant. Because the tissue is healed at the time of impression, the margin is designed to have the appropriate placement and contour. The result is a patient specific, anatomically designed, final BellaTek Abutment. Available in titanium with or without titanium nitride coating.



BellaTek Bars & Frameworks

BellaTek Bars and Frameworks are created with the technology to achieve a truly passive fit. These are available in Hader and DOLDER® Designs for overdenture bars, primary bars and in fixed hybrid designs. This technology does not require any wax or resin design from the technician. Because these are milled from one solid piece of lightweight titanium, the potential for weakness caused by soldering or laser welded joints is eliminated. You receive a completed superstructure ready for prosthesis processing. No polishing or retention creation necessary. **COMET 3i** can even tap holes for attachments.

Copymilled Bars are also available for those who prefer to design a bar or framework in a resin pattern to the patient's specifications.

NOTE: Please ensure that as many implants as necessary are used for a fully stable restoration.

Introduction



Bars & Frameworks Design Options


BellaTek® Bars and Frameworks provide simplified laboratory procedures for implant overdentures.

Hybrid Bar Design #1




- Fixed prosthesis
- Metal lingual
- Low buccal finish line
- Titanium tissue surface

Hybrid Bar Design #2



- Fixed prosthesis
- Metal lingual
- Low buccal finish line
- Titanium tissue surface

Wraparound Bar



- Fixed prosthesis
- Acrylic resin wraparound design
- Adjustments possible to accommodate for changing tissue contours
- Easy to repair and modify the acrylic resin portion of the prosthesis

Freeform Bar



- Fixed prosthesis
- Acrylic resin wraparound design

Canada Bar



- Fixed hybrid bar
- Low lingual finish line
- Titanium tissue surface
- Acrylic resin is processed to the bar

DOLDER® Egg Shape Bar



- Removable overdenture design
- Anterior clips allow for anterior/posterior rotation of the overdenture
- Provides lateral stability to the overdenture

Introduction



Bars & Frameworks Design Options (Cont'd)

DOLDER® U Shape Bar



- Removable overdenture design
- Anterior clips allow for anterior/posterior rotation of the overdenture
- Provides lateral stability to the overdenture

Hader Bar



- Removable overdenture design with clips that engage undercuts
- Provides lateral stability to the overdenture
- 4 mm between implants is needed for clips

Primary Bar




- Most secure removable overdenture design
- May be used with attachments or secondary casting
- Bar can be designed with a 2° – 6° taper

Combination Primary Bar With Hader Or DOLDER Designs





- Combines a milled bar with Hader or DOLDER Designs
- Provides excellent lateral stability to the overdenture
- Bar can be designed with a 2° – 6° taper

Design Features



Implant Cylinder Extensions



Stability Struts

Ask for these design features on the BellaTek Bars & Frameworks Work Order Form (ART880).


Copymilled Bar/Framework



- Ideal for both fixed and removable prosthetic designs
- Technician-fabricated resin pattern for porcelain or acrylic resin bonded to titanium

Introduction

Restorative Products Flowchart



Provisional Restorations

Certain® Internal Connection


PreFormance® Posts

Indications:

- Single and multiple-unit cement-retained provisional restorations
- Minimum interarch space of 7 mm
- Maximum angulation of 15°
- Intraoral use limited to 180 days
- Immediate nonocclusal loading of single-unit provisional restorations
- Multiple-unit restorations will require an evaluation of occlusal forces in lateral and protrusive excursions and masticatory forces to minimize the load on the provisional restorations
- Guided soft-tissue healing of single and multiple-unit restorations of integrated implants

Material:

PEEK: Titanium Alloy Connection



3.4 mm Seating Surface			
Emergence Profile	Collar Height	Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
3.8 mm	4 mm	IMPFP34	IMPAPF34
3.8 mm	6 mm	IMPFP36	IMPAPF36
5 mm	4 mm	IMPFP354	IMPAPF354
5 mm	6 mm	IMPFP356	IMPAPF356
4.1 mm Seating Surface			
Emergence Profile	Collar Height	Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
5 mm	4 mm	IPFP454	IPAPF454
5 mm	6 mm	IPFP456	IPAPF456
5 mm Seating Surface			
Emergence Profile	Collar Height	Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
5 mm	4 mm	IWPFP554	IPAPF554
5 mm	6 mm	IWPFP556	IPAPF556
6 mm Seating Surface			
Emergence Profile	Collar Height	Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
6 mm	4 mm	IWPFP664	IPAPF664
6 mm	6 mm	IWPFP666	IPAPF666

All PreFormance Posts are packaged with a Titanium Hexed Screw (IUNIHT).

All indicated measurements are in millimeters.

Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Provisional Restorations

Certain® Internal Connection

PreFormance® Temporary Cylinders

Indications:


- Single and multiple-unit screw-retained provisional restorations
- Minimum interarch space of 6 mm
- Maximum angulation of 15°
- Intraoral use limited to 180 days
- Immediate nonocclusal loading of single-unit provisional restorations
- Multiple-unit restorations will require an evaluation of prosthesis support to minimize the load on the provisional abutments
- Occlusal loading of single and multiple-unit restorations of integrated implants for guided tissue healing

Material:

PEEK: Titanium Alloy Connection

Description	Item Numbers
3.4 mm Seating Surface	
Hexed	IMPFTC34
Non-Hexed	IMPFTC32
4.1 mm Seating Surface	
Hexed	IPFTC41
Non-Hexed	IPFTC42
5 mm Seating Surface	
Hexed	IWPFTC51
Non-Hexed	IWPFTC52
6 mm Seating Surface	
Hexed	IWPFTC61
Non-Hexed	IWPFTC62

All PreFormance Temporary Cylinders are packaged with a Titanium Hexed Screw.



Titanium Temporary Cylinders

Indications:


- Single and multiple-unit screw-retained provisional restorations
- Minimum interarch space of 6 mm
- Maximum angulation of 15°
- Immediate nonocclusal loading of single-unit provisional restorations
- Multiple-unit restorations will require an evaluation of prosthesis support to minimize the load on the provisional abutments
- Occlusal loading of single and multiple-unit restorations of integrated implants for guided tissue healing

Material:

Titanium Alloy Connection

Description	Item Numbers
3.4 mm Seating Surface	
Hexed	IMMTCS1
Non-Hexed	—
4.1 mm Seating Surface	
Hexed	IITCS41
Non-Hexed	IITCS42
5 mm Seating Surface	
Hexed	IWTCS51
Non-Hexed	IWTCS52
6 mm Seating Surface	
Hexed	IWTCS61
Non-Hexed	IWTCS62

All Titanium Temporary Cylinders are packaged with a Titanium Hexed Screw.



Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Cement-Retained Restorations

Certain® Internal Connection


GingiHue® Posts

Indications:

- Single and multiple-unit porcelain fused-to-metal restorations
- Areas of thin tissue where gingival discoloration is possible
- Preparable to follow gingival contours
- Conventional crown and bridge procedures
- Maximum angulation correction of 15°
- Minimum interarch space of 7 mm

Material:

Titanium Alloy



3.4 mm Seating Surface

Emergence Profile	Collar Height	Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
3.8 mm	2 mm	IMAP32G	IMPAP32G
3.8 mm	4 mm	IMAP34G	IMPAP34G
5 mm	2 mm	IMAP352G	IMPAP352G
5 mm	4 mm	IMAP354G	IMPAP354G

4.1 mm Seating Surface

Emergence Profile	Collar Height	Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
5 mm	2 mm	IAPP452G	IPAP452G
5 mm	4 mm	IAPP454G	IPAP454G
6 mm	2 mm	IAPP462G	IPAP462G
6 mm	4 mm	IAPP464G	IPAP464G

5 mm Seating Surface

Emergence Profile	Collar Height	Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
5 mm	2 mm	IWPP552G	IPAP552G
5 mm	4 mm	IWPP554G	IPAP554G
6 mm	2 mm	IWPP562G	IPAP562G
6 mm	4 mm	IWPP564G	IPAP564G

6 mm Seating Surface

Emergence Profile	Collar Height	Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
6 mm	2 mm	IWPP662G	IPAP662G
6 mm	4 mm	IWPP664G	IPAP664G






All indicated measurements are in millimeters.

Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Cement-Retained Restorations

Certain® Internal Connection

GingiHue® Post Components

Pick-Up Impression Coping	Twist Lock™ Impression Coping	Laboratory Analog	Abutment Screw	Laboratory Abutment Holder
				

● = Requires the use of Hex Driver PHD02N, PHD03N or Driver Tip RASH3N, RASH8N

○ = Requires the use of Impression Coping Driver ICD00

Description	3.4 mm			4.1 mm			5 mm			6 mm	
Emergence Profile	3.4	3.8	5.0	4.1	5.0	6.0	5.0	5.6	6.0	6.0	6.8
Pick-Up Impression Coping	IMIC30	IMIC33	IMIC35	IIIC41 IIIC42*	IIIC12	IIIC60	IWIP50	IWIP55	IWIP56	IWIP60	IWIP66
Twist Lock Impression Coping		IMIT33	IMIT35	IIIC44	IIIC45	IIIC46	IWIT50	IWIT55	IWIT56	IWIT60	IWIT66
Laboratory Analog	IMMILA			IILA20			IILAW5			IILAW6	
Abutment Screw	IUNIHG										
Hexed Try-In Screw (5-pack)	IUNITS										
Laboratory Abutment Holder	ILTAH57										

*Non-Hexed

Some products may not be available outside of the United States.
Algunos productos no se comercializan fuera de EE.UU.
Certains produits ne sont pas disponibles en dehors des États-Unis.
Einige Produkte sind außerhalb der USA nicht erhältlich.

Cement-Retained Restorations

Certain® Internal Connection


UCLA Abutments


Indications:

- Laboratory fabricated custom abutments
- Single and multiple-unit cement-retained restorations
- Minimum interarch space of 6 mm
- For aesthetic restorations when tissue is limited
- Maximum angulation correction of 30°

Material:

Machined Gold Alloy Cylinder With Plastic Unitube

	
Description	Item Numbers
3.4 mm Seating Surface	
Gold (Hexed)	IMUCG1C
4.1 mm Seating Surface	
Gold (Hexed)	IGUCA1C
5 mm Seating Surface	
Gold (Hexed)	IWGA51C
6 mm Seating Surface	
Gold (Hexed)	IWGA61C








All indicated measurements are in millimeters.

Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Cement-Retained Restorations

Certain® Internal Connection

UCLA Abutment Components

Pick-Up Impression Coping	Twist Lock™ Impression Coping	Laboratory Analog	Abutment Screw	Laboratory Abutment Holder
				

● = Requires the use of Hex Driver PHD02N, PHD03N or Driver Tip RASH3N, RASH8N

○ = Requires the use of Impression Coping Driver ICD00

Description	3.4 mm			4.1 mm			5 mm			6 mm	
Emergence Profile	3.4	3.8	5.0	4.1	5.0	6.0	5.0	5.6	6.0	6.0	6.8
Pick-Up Impression Coping	IMIC30	IMIC33	IMIC35	IIIC41 IIIC42*	IIIC12	IIIC60	IWIP50	IWIP55	IWIP56	IWIP60	IWIP66
Twist Lock Impression Coping		IMIT33	IMIT35	IIIC44	IIIC45	IIIC46	IWIT50	IWIT55	IWIT56	IWIT60	IWIT66
Laboratory Analog		IMMILA		IILA20			IILAW5			IILAW6	
Abutment Screw	IUNIHG										
Hexed Try-In Screw (5-pack)	IUNITS										
Laboratory Abutment Holder	ILTAH57										

Gold Alloy Cylinder Formulation	
Melting Range	1400–1490°C (2550–2710°F)
Solidus	1400°C
Liquidus	1490°C
CTE	12.2 x 10 ⁻⁶ °K at 500°C

Gold Alloy Cylinder Melting Temperatures		
	Fahrenheit	Celsius
Melting Range	2552-2732°F	1400-1500°C
Solidus	2552-2579°F	1400-1415°C
Liquidus	2714-2732°F	1490-1500°C
CTE	12.2 x 10 ⁻⁶ °K at 932°F	12.2 x 10 ⁻⁶ °K at 500°C

*Non-Hexed

Some products may not be available outside of the United States.
Algunos productos no se comercializan fuera de EE.UU.
Certains produits ne sont pas disponibles en dehors des États-Unis.
Einige Produkte sind außerhalb der USA nicht erhältlich.

Screw-Retained Restorations

Certain® Internal Connection

UCLA Abutments


Indications:

- Single and multiple-unit restorations
- Minimum interarch space of 4 mm
- For aesthetic restorations when tissue is limited
- Non-hexed abutments packaged with screw

Material:

Machined Gold Alloy Cylinder With Plastic Unitube

Description	Item Numbers
3.4 mm Seating Surface	
Gold (Hexed)	IMUCG1C*
Gold (Non-Hexed)	IMUCG2C
Gold (Non-Hexed)	IMUCG2T
Castable (Non-Hexed)	IMUCC2C**
Castable (Non-Hexed)	IMUCC2T**
4.1 mm Seating Surface	
Gold (Hexed)	IGUCA1C*
Gold (Non-Hexed)	IGUCA2C
Gold (Non-Hexed)	IGUCA2T
Gold Standard ZR (Hexed)	—
Castable (Hexed)	—
Castable (Non-Hexed)	IUNAB2C**
Castable (Non-Hexed)	IUNAB2T**
Castable (Hexed) 25-pack	—
Castable (Non-Hexed) 25-pack	—
5 mm Seating Surface	
Gold (Hexed)	IWGA51C*
Gold (Non-Hexed)	IWGA52C
Gold (Non-Hexed)	IWGA52T
Gold Standard ZR (Hexed)	—
Castable (Hexed)	—
Castable (Non-Hexed)	IWPC52C**
Castable (Non-Hexed)	IWPC52T**
6 mm Seating Surface	
Gold (Hexed)	IWGA61C*
Gold (Non-Hexed)	IWGA62C
Gold (Non-Hexed)	IWGA62T
Gold Standard ZR (Hexed)	—
Castable (Hexed)	—
Castable (Non-Hexed)	IWPC62C**
Castable (Non-Hexed)	IWPC62T**



All indicated measurements are in millimeters.

*UCLA Gold Hexed Abutments do not include the screw. The final Gold-Tite® Screw (IUNIHG) or the Titanium Screw (IUNIHT) must be ordered separately.






**Not available in the US.

Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Screw-Retained Restorations

Certain® Internal Connection

UCLA Abutment Components

Pick-Up Impression Coping	Twist Lock™ Impression Coping	Laboratory Analog	Abutment Screw	Laboratory Abutment Holder
				

● = Requires the use of Hex Driver PHD02N, PHD03N or Driver Tip RASH3N, RASH8N

○ = Requires the use of Impression Coping Driver ICD00

Description	3.4 mm			4.1 mm			5 mm			6 mm	
Emergence Profile	3.4	3.8	5.0	4.1	5.0	6.0	5.0	5.6	6.0	6.0	6.8
Pick-Up Impression Coping	IMIC30	IMIC33	IMIC35	IIIC41 IIIC42*	IIIC12	IIIC60	IWIP50	IWIP55	IWIP56	IWIP60	IWIP66
Twist Lock Impression Coping		IMIT33	IMIT35	IIIC44	IIIC45	IIIC46	IWIT50	IWIT55	IWIT56	IWIT60	IWIT66
Laboratory Analog	IMMILA			IILA20			IILAW5			IILAW6	
Abutment Screw	IUNIHG, ILRGHG										
Hexed Try-In Screw (5-pack)	IUNITS										
Laboratory Abutment Holder	ILTAH57										

Gold Alloy Cylinder Formulation	
Melting Range	1400–1490°C (2550–2710°F)
Solidus	1400°C
Liquidus	1490°C
CTE	12.2 x 10 ⁻⁶ /°K at 500°C

*Non-Hexed

Some products may not be available outside of the United States.
Algunos productos no se comercializan fuera de EE.UU.
Certains produits ne sont pas disponibles en dehors des États-Unis.
Einige Produkte sind außerhalb der USA nicht erhältlich.

Screw-Retained Restorations

Certain® Internal Connection


Low Profile Abutments

Indications:

- Single and multiple-unit screw-retained restorations
- Minimum interarch space of 7.5 mm
- Minimum tissue height of 1 mm
- Angle correction up to 30°

Material:

Titanium Alloy



3.4 mm Seating Surface		Item Numbers			
Collar Height		Straight (Non-Hexed)	Straight (Hexed)	17° Pre-Angled (Hexed)	30° Pre-Angled (Hexed)
1 mm		ILPC341U	ILPC341	—	—
2 mm		ILPC342U	ILPC342	ILPAC3217	—
3 mm		ILPC343U	ILPC343	—	ILPAC3330
4 mm		ILPC344U	ILPC344	ILPAC3417	—
5 mm		—	—	—	ILPAC3530
4.1 mm Seating Surface		Item Numbers			
Collar Height		Straight (Non-Hexed)	Straight (Hexed)	17° Pre-Angled (Hexed)	30° Pre-Angled (Hexed)
1 mm		ILPC441U	ILPC441	—	—
2 mm		ILPC442U	ILPC442	ILPAC4217	—
3 mm		ILPC443U	ILPC443	—	ILPAC4330
4 mm		ILPC444U	ILPC444	ILPAC4417	—
5 mm		—	—	—	ILPAC4530
5 mm Seating Surface		Item Numbers			
Collar Height		Straight (Non-Hexed)	Straight (Hexed)	17° Pre-Angled (Hexed)	30° Pre-Angled (Hexed)
1 mm		ILPC541U	ILPC541	—	—
2 mm		ILPC542U	ILPC542	ILPAC5217	—
3 mm		ILPC543U	ILPC543	—	ILPAC5330
4 mm		ILPC544U	ILPC544	ILPAC5417	—
5 mm		—	—	—	ILPAC5530











All indicated measurements are in millimeters.

Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Screw-Retained Restorations

Certain® Internal Connection

Low Profile Abutment Components

<i>Pick-Up Impression Coping</i>	<i>Twist Lock™ Impression Coping</i>	<i>Laboratory Analog</i>	<i>Healing Cap</i>	<i>Gold And Castable Cylinder</i>
				
<i>QuickBridge®</i>	<i>Temporary Cylinder</i>	<i>Polishing Protector</i>	<i>Retaining Screw</i>	<i>Driver And Driver Tip</i>
				

● = Requires the use of Hex Driver PHD02N, PHD03N or Driver Tip RASH3N, RASH8N

○ = Requires the use of Impression Coping Driver ICD00

Description		4.8 mm Restorative Platform
Pick-Up Impression Coping	hexed ➤	LPCPIC1
	non-hexed ➤	LPCPIC2
Twist Lock Impression Coping	hexed ➤	LPCTIC1
	non-hexed ➤	LPCTIC2
QuickBridge Impression Coping		LPCRIC
Laboratory Analog		LPCLA
Temporary Cylinder	hexed ➤	LPCTC1
	non-hexed ➤	LPCTC2
Healing Cap		LPCHC
PreFormance® Temp. Cylinder	hexed ➤	LPCPTC1
	non-hexed ➤	LPCPTC2
QuickBridge		LPCQB
QuickBridge Replacement Cap		LPCQBCAP
Distal Extension		LPCDE
Gold Cylinder	hexed ➤	LPCGC1
	non-hexed ➤	LPCGC2
Castable Cylinder	hexed ➤	LPCCC1
	non-hexed ➤	LPCCC2
Laser Welded Cylinder		LPCTUCA
Waxing Screw		LPCWS
Retaining Screw		LPCGSH, LPCTSH
Lapping Tool		LPCAMI
Polishing Protector		LPCPP

Abutment Compatibility

Component	Hexed Abutment	Non-Hexed Abutment
Hexed	Yes	No
Non-Hexed	Yes	Yes

Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Overdenture Restorations

Certain® Internal Connection

LOCATOR® Abutments

Indications:

- Tissue supported removable overdentures on 2-4 implants
- Partially edentulous overdentures with 1 or more implants
- Limited interarch distance
- Angle correction with up to 40° between divergent implants

Material:

Abutment: Titanium Alloy With Gold-Colored Titanium Nitride Coating

Housing: Titanium Alloy

Males: Nylon

3.4 mm Seating Surface

Collar Height	Item Numbers
1 mm	IMLOA001
2 mm	IMLOA002
3 mm	IMLOA003
4 mm	IMLOA004
5 mm	IMLOA005
6 mm	IMLOA006



4.1 mm Seating Surface

Collar Height	Item Numbers
1 mm	ILOA001
2 mm	ILOA002
3 mm	ILOA003
4 mm	ILOA004
5 mm	ILOA005
6 mm	ILOA006



Included With Each Abutment:



Housing and Black Processing Male Final Male Light Retention Male Extra Light Retention Male White Spacer Ring





All indicated measurements are in millimeters.

Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Overdenture Restorations

Certain® Internal Connection

LOCATOR® Abutment Components

Impression Coping	Laboratory Analog	Retention Males	Core Tool / Driver And Driver Tip
			

Description	Item Numbers
Impression Coping	LAIC1
Laboratory Analog	LALA1
Extra Light Retention Males (4-pack) - Blue, 1.5 lbs/680 g	LAELM
Extended Range Males (4-pack) - Green, 4 lbs/1800 g	LAERM
Extra Light Retention Extended Range Males - Red, 1.5 lbs/680 g	LELARM
Replacement Males (4-pack) - White, 5 lbs/2250 g	LARMS
Light Retention Replacement Males (4-pack) - Pink, 3lbs/1350 g	LLRMS
Replacement Housing	LOAH
Replacement Housing Kit	LORHK*
Core Tool/Abutment Driver	LCTDR1
Driver Tip 24 mm(L)	LOADT4
Driver Tip 30 mm(L)	LOADT9
1.25 mm/.050 in. Driver Tip 24 mm(L) - used with LCTDR1	RASH4
1.25 mm/.050 in. Driver Tip 30 mm(L) - used with LCTDR1	RASH9
LOCATOR Bar Attachment	LOAB
LOCATOR Bar Attachment Tap	LOAT
LOCATOR Bar Attachment Drill	LOAD

*Included With LORHK
Replacement Housing Kit:



Some products may not be available outside of the United States.
Algunos productos no se comercializan fuera de EE.UU.
Certains produits ne sont pas disponibles en dehors des États-Unis.
Einige Produkte sind außerhalb der USA nicht erhältlich.

Provisional Restorations

External Hex Connection

PreFormance® Posts

Indications:

- Single and multiple-unit provisional cement-retained provisional restorations
- Minimum interarch space of 7 mm
- Maximum angulation of 15°
- Intraoral use limited to 180 days
- Immediate nonocclusal loading of single-unit provisional restorations
- Multiple-unit restorations will require an evaluation of occlusal forces in lateral and protrusive excursions and masticatory forces to minimize the load on the provisional restorations
- Guided soft-tissue healing of single and multiple-unit restorations of integrated implants

Material:

PEEK: Titanium Alloy Connection



3.4 mm Seating Surface		Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
Emergence Profile	Collar Height		
3.8 mm	4 mm	MPFP34	MPAPF34
3.8 mm	6 mm	MPFP36	MPAPF36
4.1 mm Seating Surface		Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
Emergence Profile	Collar Height		
5 mm	4 mm	PFP454	PAPF454
5 mm	6 mm	PFP456	PAPF456
5 mm Seating Surface		Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
Emergence Profile	Collar Height		
5 mm	4 mm	WPFP554	PAPF554
5 mm	6 mm	WPFP556	PAPF556
6 mm Seating Surface		Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
Emergence Profile	Collar Height		
6 mm	4 mm	WPFP664	PAPF664
6 mm	6 mm	WPFP666	PAPF666

All PreFormance Posts are packaged with a Titanium Hexed Screw (UNIHT).

All indicated measurements are in millimeters.

Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Provisional Restorations

External Hex Connection

PreFormance® Temporary Cylinders

Indications:

- Single and multiple-unit provisional screw-retained restorations
- Minimum interarch space of 6 mm
- Maximum angulation of 15°
- Intraoral use limited to 180 days
- Immediate nonocclusal loading of single-unit provisional restorations
- Multiple-unit restorations will require an evaluation of prosthesis support to minimize the load on the provisional abutments
- Occlusal loading of single and multiple-unit restorations of integrated implants for guided tissue healing

Material:

PEEK: Titanium Alloy Connection

Description	Item Numbers
3.4 mm Seating Surface	
Hexed	MPFTC34
Non-Hexed	MPFTC32
4.1 mm Seating Surface	
Hexed	PFTC41
Non-Hexed	PFTC42
5 mm Seating Surface	
Hexed	WPFTC51
Non-Hexed	WPFTC52
6 mm Seating Surface	
Hexed	WPFTC61
Non-Hexed	WPFTC62



All PreFormance Cylinders are packaged with a Titanium Hexed Screw.

Titanium Temporary Cylinders

Indications:

- Single and multiple-unit provisional screw-retained restorations
- Minimum interarch space of 6 mm
- Maximum angulation of 15°
- Immediate nonocclusal loading of single-unit provisional restorations
- Multiple-unit restorations will require an evaluation of prosthesis support to minimize the load on the provisional abutments
- Occlusal loading of single and multiple-unit restorations of integrated implants for guided tissue healing

Material:

Titanium Alloy Connection

Description	Item Numbers
3.4 mm Seating Surface	
Hexed	MMTCS1
Non-Hexed	—
4.1 mm Seating Surface	
Hexed	ITCS41
Non-Hexed	ITCS42
5 mm Seating Surface	
Hexed	WTCS51
Non-Hexed	WTCS52
6 mm Seating Surface	
Hexed	WTCS61
Non-Hexed	WTCS62



All Titanium Temporary Cylinders are packaged with a Titanium Screw. (IUNIHT).

All indicated measurements are in millimeters.

Some products may not be available outside of the United States.

Cement-Retained Restorations

External Hex Connection



GingiHue® Posts

Indications:

- Single and multiple-unit porcelain fused-to-metal restorations
- Areas of thin tissue where gingival discoloration is possible
- Preparable to follow gingival contours
- Conventional crown and bridge procedures
- Maximum angulation correction of 15°
- Minimum interarch space of 7 mm

Material:

Straight Post (Commercially Pure Titanium)
Pre-Angled Post (Titanium Alloy)



3.4 mm Seating Surface			
Emergence Profile	Collar Height	Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
3.8 mm	2 mm	MAP32G	MPAP32G
3.8 mm	4 mm	MAP34G	MPAP34G
4.1 mm Seating Surface			
Emergence Profile	Collar Height	Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
5 mm	2 mm	APP452G	PAP452G
5 mm	4 mm	APP454G	PAP454G
6 mm	2 mm	APP462G	PAP462G
6 mm	4 mm	APP464G	PAP464G
5 mm Seating Surface			
Emergence Profile	Collar Height	Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
5 mm	2 mm	WPP552G	PAP552G
5 mm	4 mm	WPP554G	PAP554G
6 mm	2 mm	WPP562G	PAP562G
6 mm	4 mm	WPP564G	PAP564G
6 mm Seating Surface			
Emergence Profile	Collar Height	Item Numbers (Straight)	Item Numbers (15° Pre-Angled)
6 mm	2 mm	WPP662G	PAP662G
6 mm	4 mm	WPP664G	PAP664G










All indicated measurements are in millimeters.

Some products may not be available outside of the United States.
Algunos productos no se comercializan fuera de EE.UU.
Certains produits ne sont pas disponibles en dehors des États-Unis.
Einige Produkte sind außerhalb der USA nicht erhältlich.

Cement-Retained Restorations

External Hex Connection

GingiHue® Post Components

Pick-Up Impression Coping	Twist Lock™ Impression Coping	Laboratory Analog	Abutment Screw	Laboratory Abutment Holder
				
				

● = Requires the use of Hex Driver PHD02N, PHD03N or Driver Tip RASH3N, RASH8N

○ = Requires the use of Impression Coping Driver ICD00

■ = Requires the use of Square Driver PSQD0N, PSQD1N or Driver Tip RASQ3N, RASQ8N

Description	3.4 mm	4.1 mm	5 mm	6 mm
<i>Emergence Profile</i>	3.8	5.0 6.0	5.0 6.0	6.0
Pick-Up Impression Coping	MIC33	IIC12 IIC60	WIP55 WIP56	WIP66
Twist Lock Impression Coping	—	IIC45 IIC46	WIT55 WIT56	WIT66
Laboratory Analog	MMILA	ILA20	ILAW5	ILAW6
Abutment Screw	UNISG*, UNIHG			
Square Try-In Screw (5-pack)	MUNITS	UNITS		
Laboratory Abutment Holder	LTAH5	LTAH7		

* Recommended

Some products may not be available outside of the United States.
Algunos productos no se comercializan fuera de EE.UU.
Certains produits ne sont pas disponibles en dehors des États-Unis.
Einige Produkte sind außerhalb der USA nicht erhältlich.

Cement-Retained Restorations

External Hex Connection


UCLA Abutments

Indications:

- Laboratory fabricated custom abutments
- Single and multiple-unit cement retained restorations
- Minimum interarch space of 6 mm
- For aesthetic restorations when tissue is limited
- Maximum angulation correction of 30°

Material:

Machined Gold Alloy Cylinder With Plastic Unitube

	
Description	Item Numbers
3.4 mm Seating Surface	
Gold (Hexed)	MUCG1C
4.1 mm Seating Surface	
Gold (Hexed)	GUCA1C
Gold Standard ZR (Hexed)	SGUCA1C
Castable (Hexed)	UNAB1C*
Castable (Hexed) 25-pack	UNA125*
5 mm Seating Surface	
Gold (Hexed)	WGA51C
Gold Standard ZR (Hexed)	SWGA51C
Castable (Hexed)	WPC51C*
6 mm Seating Surface	
Gold (Hexed)	WGA61C
Gold Standard ZR (Hexed)	SWGA61C
Castable (Hexed)	WPC61C*



All indicated measurements are in millimeters.










*Not available in the US.

Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Cement-Retained Restorations

External Hex Connection

UCLA Abutment Components

Pick-Up Impression Coping	Twist Lock™ Impression Coping	Laboratory Analog	Abutment Screw	Laboratory Abutment Holder
				
				
● = Requires the use of Hex Driver PHD02N, PHD03N or Driver Tip RASH3N, RASH8N		○ = Requires the use of Impression Coping Driver ICD00	■ = Requires the use of Square Driver PSQD0N, PSQD1N or Driver Tip RASQ3N, RASQ8N	

Description	3.4 mm	4.1 mm		5 mm		6 mm
Emergence Profile	3.8	5.0	6.0	5.0	6.0	6.0
Pick-Up Impression Coping	MIC33	IIC12	IIC60	WIP55	WIP56	WIP66
Twist Lock Impression Coping	—	IIC45	IIC46	WIT55	WIT56	WIT66
Laboratory Analog	MMILA	ILA20		ILAW5		ILAW6
Abutment Screw	UNISG*, UNIHG					
Square Try-In Screw (5-pack)	MUNITS	UNITS				
Laboratory Abutment Holder	LTAH5	LTAH7				

Gold Alloy Cylinder Formulation	
Melting Range	1400–1490°C (2550–2710°F)
Solidus	1400°C
Liquidus	1490°C
CTE	12.2 x 10 ⁻⁶ °K at 500°C

Gold Alloy Cylinder Melting Temperatures		
	Fahrenheit	Celsius
Melting Range	2552-2732°F	1400-1500°C
Solidus	2552-2579°F	1400-1415°C
Liquidus	2714-2732°F	1490-1500°C
CTE	12.2 x 10 ⁻⁶ °K at 932°F	12.2 x 10 ⁻⁶ °K at 500°C

*Recommended

Some products may not be available outside of the United States.
Algunos productos no se comercializan fuera de EE.UU.
Certains produits ne sont pas disponibles en dehors des États-Unis.
Einige Produkte sind außerhalb der USA nicht erhältlich.

Screw-Retained Restorations

External Hex Connection


UCLA Abutments


Indications:

- Single and multiple-unit restorations
- Minimum interarch space of 4 mm
- For aesthetic restorations when tissue is limited

Material:

Machined Gold Alloy Cylinder With Plastic Unitube

	
Description	Item Numbers
3.4 mm Seating Surface	
Gold (Hexed)	MUCG1C
Gold (Non-Hexed)	MUCG2C
Gold (Non-Hexed)	—
Castable (Hexed)	MUCC1C*
Castable (Non-Hexed)	MUCC2C*
Castable (Non-Hexed)	—
4.1 mm Seating Surface	
Gold (Hexed)	GUCA1C
Gold (Non-Hexed)	GUCA2C
Gold (Non-Hexed)	—
Gold Standard ZR (Hexed)	SGUCA1C
Castable (Hexed)	UNAB1C*
Castable (Non-Hexed)	UNAB2C*
Castable (Non-Hexed)	—
Castable (Hexed) 25-pack	UNA125*
Castable (Non-Hexed) 25-pack	UNA225*
5 mm Seating Surface	
Gold (Hexed)	WGA51C
Gold (Non-Hexed)	WGA52C
Gold (Non-Hexed)	—
Gold Standard ZR (Hexed)	SWGA51C
Castable (Hexed)	WPC51C*
Castable (Non-Hexed)	WPC52C*
Castable (Non-Hexed)	—
6 mm Seating Surface	
Gold (Hexed)	WGA61C
Gold (Non-Hexed)	WGA62C
Gold (Non-Hexed)	—
Gold Standard ZR (Hexed)	SWGA61C
Castable (Hexed)	WPC61C*
Castable (Non-Hexed)	WPC62C*
Castable (Non-Hexed)	—



All indicated measurements are in millimeters.










*Not available in the US.

Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Screw-Retained Restorations

External Hex Connection

UCLA Abutment Components

Pick-Up Impression Coping	Twist Lock™ Impression Coping	Laboratory Analog	Abutment Screw	Laboratory Abutment Holder
				
				

● = Requires the use of Hex Driver PHD02N, PHD03N or Driver Tip RASH3N, RASH8N

○ = Requires the use of Impression Coping Driver ICD00

■ = Requires the use of Square Driver PSQD0N, PSQD1N or Driver Tip RASQ3N, RASQ8N

Description	3.4 mm	4.1 mm		5 mm		6 mm
Emergence Profile	3.8	5.0	6.0	5.0	6.0	6.0
Pick-Up Impression Coping	MIC33	IIC12	IIC60	WIP55	WIP56	WIP66
Twist Lock Impression Coping	—	IIC45	IIC46	WIT55	WIT56	WIT66
Laboratory Analog	MMILA	ILA20		ILAW5		ILAW6
Abutment Screw	UNISG*, UNIHG					
Square Try-In Screw (5-pack)	MUNITS	UNITS				
Laboratory Abutment Holder	LTAH5	LTAH7				

Gold Alloy Cylinder Melting Temperatures		
	Fahrenheit	Celsius
Melting Range	2552-2732°F	1400-1500°C
Solidus	2552-2579°F	1400-1415°C
Liquidus	2714-2732°F	1490-1500°C
CTE	12.2 x 10 ⁻⁶ °K at 932°F	12.2 x 10 ⁻⁶ °K at 500°C

* Recommended

Some products may not be available outside of the United States.
Algunos productos no se comercializan fuera de EE.UU.
Certains produits ne sont pas disponibles en dehors des États-Unis.
Einige Produkte sind außerhalb der USA nicht erhältlich.

Screw-Retained Restorations

External Hex Connection


Low Profile Abutments

Indications:

- Single and multiple-unit screw-retained restorations
- Minimum interarch space of 7.5 mm
- Minimum tissue height of 1 mm
- Angle correction up to 30°
- External hex connection 3.4 mm(D) Low Profile Abutments are limited for use in the anterior

Material:

Titanium Alloy

				
3.4 mm Seating Surface				
Collar Height	Item Numbers			
	Straight (Non-Hexed)	Straight (Hexed)	17° Pre-Angled (Hexed)	30° Pre-Angled (Hexed)
1 mm	LPC341U	LPC341	—	—
2 mm	LPC342U	LPC342	LPAC3217	—
3 mm	LPC343U	LPC343	—	LPAC3330
4 mm	LPC344U	LPC344	LPAC3417	—
5 mm	—	—	—	LPAC3530
4.1 mm Seating Surface				
Collar Height	Item Numbers			
	Straight (Non-Hexed)	Straight (Hexed)	17° Pre-Angled (Hexed)	30° Pre-Angled (Hexed)
1 mm	LPC441U	LPC441	—	—
2 mm	LPC442U	LPC442	LPAC4217	—
3 mm	LPC443U	LPC443	—	LPAC4330
4 mm	LPC444U	LPC444	LPAC4417	—
5 mm	—	—	—	LPAC4530
5 mm Seating Surface				
Collar Height	Item Numbers			
	Straight (Non-Hexed)	Straight (Hexed)	17° Pre-Angled (Hexed)	30° Pre-Angled (Hexed)
1 mm	LPC541U	LPC541	—	—
2 mm	LPC542U	LPC542	LPAC5217	—
3 mm	LPC543U	LPC543	—	LPAC5330
4 mm	LPC544U	LPC544	LPAC5417	—
5 mm	—	—	—	LPAC5530











All indicated measurements are in millimeters.

Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Screw-Retained Restorations

External Hex Connection

Low Profile Abutment Components

<i>Pick-Up Impression Coping</i>	<i>Twist Lock™ Impression Coping</i>	<i>Laboratory Analog</i>	<i>Healing Cap</i>	<i>Gold And Castable Cylinder</i>
				
<i>QuickBridge®</i>	<i>Temporary Cylinder</i>	<i>Polishing Protector</i>	<i>Retaining Screw</i>	<i>Driver And Driver Tip</i>
				

● = Requires the use of Hex Driver PHD02N, PHD03N or Driver Tip RASH3N, RASH8N

○ = Requires the use of Impression Coping Driver ICD00

Description		4.8 mm Restorative Platform
Pick-Up Impression Coping	hexed ➤	LPCPIC1
	non-hexed ➤	LPCPIC2
Twist Lock Impression Coping	hexed ➤	LPCTIC1
	non-hexed ➤	LPCTIC2
QuickBridge Impression Coping		LPCRIC
Laboratory Analog		LPCLA
Temporary Cylinder	hexed ➤	LPCTC1
	non-hexed ➤	LPCTC2
Healing Cap		LPCHC
PreFormance® Temp. Cylinder	hexed ➤	LPCPTC1
	non-hexed ➤	LPCPTC2
QuickBridge		LPCQB
QuickBridge Replacement Cap		LPCQBCAP
Distal Extension		LPCDE
Gold Cylinder	hexed ➤	LPCGC1
	non-hexed ➤	LPCGC2
Castable Cylinder	hexed ➤	LPCCC1
	non-hexed ➤	LPCCC2
Laser Welded Cylinder		LPCTUCA
Waxing Screw		LPCWS
Retaining Screw		LPCGSH, LPCTSH
Lapping Tool		LPCAMI
Polishing Protector		LPCPP

Abutment Compatibility

Component	Hexed Abutment	Non-Hexed Abutment
Hexed	Yes	No
Non-Hexed	Yes	Yes

Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Overdenture Restorations

External Hex Connection

LOCATOR® Abutments

Indications:

- Tissue supported removable overdentures on 2-4 implants
- Partially edentulous overdentures with 1 or more implants
- Limited interarch distance
- Angle correction with up to 40° between divergent implants

Material:


Abutment: Titanium Alloy With Gold-Colored Titanium Nitride Coating

Housing: Titanium Alloy

Males: Nylon


3.4 mm Seating Surface

Collar Height	Item Numbers
1 mm	MLOA001
2 mm	MLOA002
3 mm	MLOA003
4 mm	MLOA004
5 mm	MLOA005
6 mm	MLOA006




4.1 mm Seating Surface

Collar Height	Item Numbers
1 mm	LOA001
2 mm	LOA002
3 mm	LOA003
4 mm	LOA004
5 mm	LOA005
6 mm	LOA006




Included With Each Abutment:



Temporary Healing Retention Cylinder

Retains denture with soft-tissue liner.







All indicated measurements are in millimeters.

Some products may not be available outside of the United States.
Algunos productos no se comercializan fuera de EE.UU.
Certains produits ne sont pas disponibles en dehors des États-Unis.
Einige Produkte sind außerhalb der USA nicht erhältlich.

Overdenture Restorations

External Hex Connection

LOCATOR® Abutment Components

Impression Coping	Laboratory Analog	Retention Males	Core Tool / Driver And Driver Tip
			

Description	Item Numbers
Impression Coping	LAIC1
Laboratory Analog	LALA1
Extra Light Retention Males (4-pack) - Blue, 1.5 lbs/680 g	LAELM
Extended Range Males (4-pack) - Green, 4 lbs/1800 g	LAERM
Extra Light Retention Extended Range Males - Red, 1.5 lbs/680 g	LELARM
Replacement Males (4-pack) - White, 5 lbs/2250 g	LARMS
Light Retention Replacement Males (4-pack) - Pink, 3 lbs/1350 g	LLRMS
Replacement Housing	LOAH
Replacement Housing Kit	LORHK*
Core Tool/Abutment Driver	LCTDR1
Driver Tip 24 mm(L)	LOADT4
Driver Tip 30 mm(L)	LOADT9
1.25 mm/.050 in. Driver Tip 24 mm(L) - used with LCTDR1	RASH4
1.25 mm/.050 in. Driver Tip 30 mm(L) - used with LCTDR1	RASH9
LOCATOR Bar Attachment	LOAB
LOCATOR Bar Attachment Tap	LOAT
LOCATOR Bar Attachment Drill	LOAD

*Included With LORHK
Replacement Housing Kit:




Some products may not be available outside of the United States.
Algunos productos no se comercializan fuera de EE.UU.
Certains produits ne sont pas disponibles en dehors des États-Unis.
Einige Produkte sind außerhalb der USA nicht erhältlich.

Screws


Certain® Internal Connection

Abutment Screws




	<i>Gold-Tite® Screw</i>	<i>Titanium Screw</i>	<i>Large Gold-Tite Screw</i>	<i>Large Titanium Screw</i>
Screws	IUNIHG	IUNIHT	ILRGHG	ILRGHT
Drivers	PHD02N, PHD03N	PHD02N, PHD03N	PHD02N, PHD03N	PHD02N, PHD03N
Driver Tips	RASH3N, RASH8N	RASH3N, RASH8N	RASH3N, RASH8N	RASH3N, RASH8N
Torque	20 Ncm	20 Ncm	20 Ncm	20 Ncm
For Use With:	UCLA Abutment, GingiHue® Post, Hexed Temporary Cylinder	UCLA Abutment, GingiHue Post, Hexed Temporary Cylinder	Non-Hexed UCLA Abutment, Non-Hexed Temporary Cylinder, BellaTek® Bars	

Retaining Screws



	<i>Low Profile Abutment Gold-Tite Screw</i>	<i>Gold-Tite Screw</i>	<i>Low Profile Abutment Titanium Screw</i>	<i>Titanium Hexed Screw</i>
Screws	LPCGSH	GSH20, GSH30, GSH70	LPCTSH	TSH30
Drivers	PHD02N, PHD03N	PHD02N, PHD03N	PHD02N, PHD03N	PHD02N, PHD03N
Driver Tips	RASH3N, RASH8N	RASH3N, RASH8N	RASH3N, RASH8N	RASH3N, RASH8N
Torque	10 Ncm	10 Ncm	10 Ncm	10 Ncm
For Use With:	Low Profile Abutment	GSH20 - Standard Abutment 6 mm(D) X 2 mm(H) Conical Abutment 6 mm(D) X 1 mm(H) GSH30 - Standard, Conical and IOL® Abutments GSH70 - Standard, Conical and IOL Abutments	Low Profile Abutment	Standard, Conical and IOL Abutments

Waxing Screws & Try-In Screws




	<i>Laboratory Waxing Screw</i>	<i>Hexed Try-In Screw (5-pack)</i>
Screws	IWSU30	IUNITS
Drivers	PHD02N, PHD03N	PHD02N, PHD03N
Driver Tips	RASH3N, RASH8N	RASH3N, RASH8N
For Use With:	Pick-Up and Impression Copings, Implant Level Laboratory Procedures	Hexed UCLA Abutment, GingiHue® Post, GingiHue Post 15°

Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Screws


External Hex Connection

Abutment Screws




	<i>Gold-Tite® Square</i>	<i>Gold-Tite Hexed</i>	<i>Titanium Hexed</i>
Screws	UNISG	UNIHG	UNIHT
Drivers	PSQD0N, PSQD1N	PHD02N, PHD03N	PHD02N, PHD03N
Driver Tips	RASQ3N, RASQ8N	RASH3N, RASH8N	RASH3N, RASH8N
Torque	32–35 Ncm	20 Ncm	20 Ncm
For Use With:	UCLA Abutment, GingiHue® Post, GingiHue Post 15°, Titanium Cylinder		

Retaining Screws



	<i>Low Profile Abutment Gold-Tite Screw</i>	<i>Gold-Tite Screw</i>	<i>Low Profile Abutment Titanium Screw</i>	<i>Titanium Hexed Screw</i>
Screws	LPCGSH	GSH20, GSH30, GSH70	LPCTSH	TSH30
Drivers	PHD02N, PHD03N	PHD02N, PHD03N	PHD02N, PHD03N	PHD02N, PHD03N
Driver Tips	RASH3N, RASH8N	RASH3N, RASH8N	RASH3N, RASH8N	RASH3N, RASH8N
Torque	10 Ncm	10 Ncm	10 Ncm	10 Ncm
For Use With:	Low Profile Abutment	GSH20 - Standard Abutment 6 mm(D) X 2 mm(H) Conical Abutment 6 mm(D) X 1 mm(H) GSH30 - Standard, Conical and IOL® Abutments GSH70 - Standard, Conical and IOL Abutments	Low Profile Abutment	Standard, Conical and IOL Abutments

Waxing Screws & Try-In Screws



	<i>Abutment Waxing Screw/Guide Pin</i>	<i>Laboratory Waxing Screw</i>	<i>Square Try-In Screw (5-pack)</i>
Screws	WSK10, WSK15	WSK30, WSU30	UNITS, MUNITS
Drivers	PHD02N, PHD03N	PHD02N, PHD03N	PSQD0N, PSQD1N
Driver Tips	RASH3N, RASH8N	RASH3N, RASH8N	RASQ3N, RASQ8N
For Use With:	Standard and Conical Abutments	Implant Level Laboratory Procedures	UCLA Abutment, GingiHue Post, GingiHue Post 15°

Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Miscellaneous Restorative Products

Mechanical Driver Kits

NCATDO *Contra Angle Torque Drivers And Accessory Kit**



NPSDK0 *Prosthetic Drivers And Accessory Kit**



NCATDOC *Contra Angle Torque Driver And Accessory Kit Certain®**



L-TIRWK *Low Torque Indicating Ratchet Wrench And Accessory Kit***



NCATDO <i>Contra Angle Torque Drivers And Accessory Kit</i>	NCATDOC <i>Contra Angle Torque Driver And Accessory Kit Certain</i>	NPSDK0 <i>Prosthetic Drivers And Accessory Kit</i>	L-TIRWK <i>Low Torque Indicating Ratchet Wrench And Accessory Kit</i>	<i>Description</i>	<i>Item Numbers</i>
•	•	•		Narrow Posterior Hex Driver, 17 mm(L)	PHD02N
•		•		Narrow Posterior Square Driver, 17 mm(L)	PSQD0N
•	•	•		Contra Angle Torque Driver Body	HTD-C*
•	•	•	•	Narrow Right Angle Hex Driver Tip, 24 mm(L)	RASH3N
•	•	•		Narrow Right Angle Hex Driver Tip, 30 mm(L)	RASH8N
•		•	•	Narrow Right Angle Square Driver Tip, 24 mm(L)	RASQ3N
•		•		Narrow Right Angle Square Driver Tip, 30 mm(L)	RASQ8N
•	•	•		Right Angle Abutment Driver Tip - Steel	RASA3
		•		Contra Angle Torque Controller, 10 Ncm	TD-10*
•	•	•		Contra Angle Torque Controller, 20 Ncm	TD-20*
•		•		Contra Angle Torque Controller, 35 Ncm	TD-35*
•	•	•		Contra Pre-action Handle	CATDH
•	•	•		Prosthetic System Driver Tray	PSDT1
		•		Narrow Impression Coping Driver	ICD00
		•		Posterior Abutment Driver, 17 mm(L)	PAD00
		•		Standard Abutment Driver, 24 mm(L)	PAD24
		•		Posterior Slotted Screw Driver, 17 mm(L)	PSD00
		•		Standard Slotted Screw Driver, 24 mm(L)	PSD01
		•		Narrow Standard Slotted Square Driver, 24 mm(L)	PSQD1N
		•		Narrow Standard Hex Driver, 24 mm(L)	PHD03N
		•		Right Angle Slotted Driver Tip, 24 mm(L)	RASD1
		•		Right Angle Slotted Driver Tip, 30 mm(L)	RASD6
			•	Low Torque Indicating Ratchet Wrench	L-TIRW**
			•	Restorative Torque Tray	RTI2035TR

* HTD-C, TD-10, TD-20, and TD-35 are manufactured by NSK America Corp/NSK dental LLC and distributed by BIOMET 3i, LLC. Sold separately.

** L-TIRW is manufactured by Elos Medtech AB and distributed by BIOMET 3i, LLC. Sold separately.

Some products may not be available outside of the United States.
Algunos productos no se comercializan fuera de EE.UU.
Certains produits ne sont pas disponibles en dehors des États-Unis.
Einige Produkte sind außerhalb der USA nicht erhältlich.

Miscellaneous Restorative Products

Prosthodontic Screwdriver

The cordless Prosthodontic Screwdriver has torque control, which enables precise tightening of retaining screws. This cordless screwdriver reduces the need for traditional mechanical rotary ratchet wrench systems.



Features Include:

- Torque control of 8-40 Ncm in 1 Ncm steps for precise tightening of retaining screws
- An automatic stop when the preset torque is reached
- Cordless operation with constant drive speed
- A small contra-angle head for better visibility of the treatment site
- A contra-angle that can be rotated 360°
- A clockwise and counter clockwise operation
- A sterilizable contra-angle

Prosthodontic Screwdriver IA-400

Consists Of:

- A handpiece drive IA-40H with keypad and rechargeable Li-Ion battery (not sterilizable)
- A contra-angle IA-80 with a hexagon chucking system
- A charging station with country specific cables available


Orderable Parts

	<i>Item Number</i>
IA-400 Prosthodontic Screw Driver	16934000
IA-400 Charging Station	05338301
IA-80 Contra-Angle Hand Piece, 80:1	16338001
IA-40H Handpiece Drive	16834000
Operation Lever	05937000
Power Cord US, CAN, J	05456100
Power Cord UK, IRL	05457700
Power Cord EU	05464400
Power Cord AUS, NZ	05467600
Rechargeable Battery, Li-Ion	05451300


Some products may not be available outside of the United States.
Algunos productos no se comercializan fuera de EE.UU.
Certains produits ne sont pas disponibles en dehors des États-Unis.
Einige Produkte sind außerhalb der USA nicht erhältlich.

Miscellaneous Restorative Products


Drivers & Driver Tips




	Small Hex		Large Hex		Square	
Length	Driver	Tip	Driver	Tip	Driver	Tip
17 mm	PHD00N	—	PHD02N	—	PSQD0N	—
24 mm	PHD01N	RASH2N	PHD03N	RASH3N	PSQD1N	RASQ3N
30 mm	—	RASH7N	—	RASH8N	—	RASQ8N
For Use With:	External Hex Cover Screws		Healing Abutments, Abutment Screws, Certain® Cover Screws and Retaining Screws		Square Gold-Tite® Abutment Screws and Square Try-in Screws	



	Abutment		LOCATOR®
Length	Driver	Tip	Tip Only
17 mm	PAD00	—	—
24 mm	PAD24	RASA3	LOADT4
30 mm	—	—	LOADT9
For Use With:	Low Profile Abutments		LOCATOR Abutments



	Impression Coping	.050" Hex
Length	Driver Only	Tip Only
17 mm	ICD00	—
24 mm	—	RASH4
30 mm	—	RASH9
For Use With:	Twist Lock™ Transfer Impression Copings	LOCATOR Abutments with LOCATOR Core Tools


MIDT Manual Interchangeable Driver Tip Handle
 

Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Miscellaneous Restorative Products

Laboratory Tools

Lapping Tool For Castable UCLA Abutment With Guide




RH600 Reamer With Handle




IQSA01 Certain QuickSeat® Activator Tool



Polishing Protector, UCLA Abutment



Laboratory Abutment Holder



Miscellaneous Items

L-TIRW Low Torque Indicating Ratchet Wrench*

L-TIRWK Low Torque Indicating Ratchet Wrench Kit*

C9980


Standard ISO 1797 Adapter

C9981


Short ISO 1797 Adapter




AG900 Prosthetic Angle Guide Kit




TMP80 Tissue Measuring Post



Certain Screw Removal Tools



Screw Removal Tools



All indicated measurements are in millimeters.

* L-TIRW and L-TIRWK are manufactured by Elos Medtech AB and distributed by BIOMET 3i, LLC.

Some products may not be available outside of the United States.
Algunos productos no se comercializan fuera de EE.UU.
Certains produits ne sont pas disponibles en dehors des États-Unis.
Einige Produkte sind außerhalb der USA nicht erhältlich.

Miscellaneous Restorative Products


Tapered Navigator® Laboratory Kit

SGTLKIT Tapered Navigator Laboratory Kit



For Certain® Tapered Implants Including:

- Certain 3.25, 4 and 5 mm Implants
- Certain Tapered PREVAIL® 4/3 and 5/4 Implants

 SGTLKIT Surgical Kit	Description	Item Numbers
•	Tapered Navigator Laboratory Tray	SGTLTRAY
•	Navigator Certain Analog Mount, 3.4 mm(D) - Short	MSGIAMS
•	Navigator Certain Analog Mount, 3.4 mm(D) - Long	MSGIAML
•	Navigator Certain Analog Mount, 4.1 mm(D) - Short	SGIAM4S
•	Navigator Certain Analog Mount, 4.1 mm(D) - Long	SGIAM4L
•	Navigator Certain Analog Mount, 5 mm(D) - Short	SGIAM5S
•	Navigator Certain Analog Mount, 5 mm(D) - Long	SGIAM5L


Parallel Walled Navigator Laboratory Kit

SGLKIT Navigator Laboratory Kit



For Certain Parallel Walled Implants Including:

- Certain 3.25, 4 and 5 mm Implants
- Certain Tapered PREVAIL 4/3 and 5/4 Implants

 SGLKIT Laboratory Kit	Description	Item Numbers
•	Navigator Laboratory Tray	SGLTRAY
•	Navigator Certain Analog Mount, 3.4 mm(D) x 1 mm(L)	MSGIAM1
•	Navigator Certain Analog Mount, 3.4 mm(D) x 2 mm(L)	MSGIAM2
•	Navigator Certain Analog Mount, 3.4 mm(D) x 3 mm(L)	MSGIAM3
•	Navigator Certain Analog Mount, 3.4 mm(D) x 4 mm(L)	MSGIAM4
•	Navigator Certain Analog Mount, 4.1 mm(D) x 1 mm(L)	SGIAM41
•	Navigator Certain Analog Mount, 4.1 mm(D) x 2 mm(L)	SGIAM42
•	Navigator Certain Analog Mount, 4.1 mm(D) x 3 mm(L)	SGIAM43
•	Navigator Certain Analog Mount, 4.1 mm(D) x 4 mm(L)	SGIAM44
•	Navigator Certain Analog Mount, 5 mm(D) x 1 mm(L)	SGIAM51
•	Navigator Certain Analog Mount, 5 mm(D) x 2 mm(L)	SGIAM52
•	Navigator Certain Analog Mount, 5 mm(D) x 3 mm(L)	SGIAM53
•	Navigator Certain Analog Mount, 5 mm(D) x 4 mm(L)	SGIAM54

Some products may not be available outside of the United States.
Algunos productos no se comercializan fuera de EE.UU.
Certains produits ne sont pas disponibles en dehors des États-Unis.
Einige Produkte sind außerhalb der USA nicht erhältlich.

Miscellaneous Restorative Products

Navigator® System Instrumentation Overview



Tapered Navigator Surgical Kit



Tapered Navigator Laboratory Kit



Parallel Walled Navigator Surgical Kit



Parallel Walled Navigator Laboratory Kit

The Navigator System For Guided Surgery

The Navigator System is designed to help clinicians transform computer-based planning into practical implant placement with an opportunity for immediate provisionalization.

The Navigator System features control of hex-orientation to fabricate and deliver provisional restorations immediately following surgery. This unique hex-orientation feature is especially beneficial when immediately loading single tooth cases, cement-retained cases and cases where angulation correction is required.

The Navigator System enables placement and restoration of the following Biomet 3i Tapered Platform Switched and Parallel Walled Implant Systems:

- Tapered/Parallel Walled Certain® 3.25, 4 and 5 mm Implants
- Certain PREVAIL® 4/3, 5/4 Implants
- Certain Tapered Platform Switched 4/3, 5/4 and 6/5 mm Implants

An integral part of the Navigator System is the Navigator Laboratory Kit, which contains easy-to-identify components needed to fabricate a master cast and provisional prosthesis prior to the day of surgery. The Navigator Laboratory Kit provides you with the opportunity to:

- Treat patients for immediate-load cases
- Be part of a collaborative approach to comprehensive patient treatment
- Embrace emerging technology with open architecture compatibility that provides freedom to choose among leading planning software solutions

Patient Education Products

<i>Description</i>	<i>Item Numbers</i>
Patient Ed Model - DIEM®2 (4-unit)	DIEM2-4UNIT
Patient Ed Model - DIEM 2 (6-unit)	DIEM2-6UNIT
Patient Ed Model - OSSEOTITE® Implant	EXHEXMOD
Patient Ed Model - LOCATOR®	PEMODLOC
Patient Ed Model - Single/Multiple Teeth	PEMODSM
Patient Ed Brochure - Bone Loss (25-pack)	ART885
Patient Ed Brochure - Implant Retained, Fixed (25-pack)	ART888
Patient Ed Brochure - Implant Retained, LOCATOR (25-pack)	ART889
Patient Ed Brochure - Implant Retained, Bar (25-pack)	ART890
Patient Ed Brochure - Immediate Teeth Replacement, DIEM (25-pack)	ART929
Patient Ed Care Card (25-pack)	ART930
Patient Ed Flip Chart	ART943
Patient Ed DVD - Multilingual Version	ART951EUR
Patient Ed Brochure - Bone Loss (Spanish Version)	ART957S



Some products may not be available outside of the United States.
 Algunos productos no se comercializan fuera de EE.UU.
 Certains produits ne sont pas disponibles en dehors des États-Unis.
 Einige Produkte sind außerhalb der USA nicht erhältlich.

Cat.#	Description	Page#	Cat.#	Description	Page#
05338301	IA-400 CHARGING STATION	28	IIC44	CERTAIN STRAIGHT TWIST LOCK COPING 4.1 mm(D) X 4.1 mm(P)	4,6,8
05451300	RECHARGEABLE BATTERY, LI-ION	28	IIC45	CERTAIN EP TWIST LOCK COPING 4.1 mm(D) X 5 mm(P)	4,6,8
05456100	POWER CORD US, CAN, J	28	IIC46	CERTAIN EP TWIST LOCK COPING 4.1 mm(D) X 6 mm(P)	4,6,8
05457700	POWER CORD UK, IRL	28	IIC60	CERTAIN EP PICK-UP COPING 4.1 mm(D) X 6 mm(P)	4,6,8
05464400	POWER CORD EU	28	IILA20	CERTAIN IMPLANT LAB ANALOG 4.1 mm(D)	4,6,8
05467600	POWER CORD AUS, NZ	28	IILAW5	CERTAIN IMPLANT LAB ANALOG 5 mm(D)	4,6,8
05937000	OPERATION LEVER	28	IILAW6	CERTAIN IMPLANT LAB ANALOG 6 mm(D)	4,6,8
16338001	IA-80 CONTRA-ANGLE HAND PIECE, 80:1	28	IITCS41	CERTAIN TEMPORARY (HEXED) CYLINDER 4.1 mm(D) W/HEXED SCREW	2
16834000	IA-40H HANDPIECE DRIVE	28	IITCS42	CERTAIN TEMPORARY (NON-HEXED) CYLINDER 4.1 mm(D) W/LARGE DIA. SCREW	2
16934000	IA-400 PROSTHODONTIC SCREW DRIVER	28	ILA20	IMPLANT LAB ANALOG 4.1 mm(D)	16,18,20
AG900	PROSTHETIC ANGLE GUIDE KIT (9-PIECE SET)	30	ILAW5	IMPLANT LAB ANALOG 5 mm(D)	16,18,20
APP452G	GINGIHUE® POST 4.1 mm(D) X 5 mm(P) X 2 mm(H)	15	ILAW6	IMPLANT LAB ANALOG 6 mm(D)	16,18,20
APP454G	GINGIHUE POST 4.1 mm(D) X 5 mm(P) X 4 mm(H)	15	ILOA001	CERTAIN LOCATOR ABUTMENT 4.1 mm(D) X 1 mm(H)	11
APP462G	GINGIHUE POST 4.1 mm(D) X 6 mm(P) X 2 mm(H)	15	ILOA002	CERTAIN LOCATOR ABUTMENT 4.1 mm(D) X 2 mm(H)	11
APP464G	GINGIHUE POST 4.1 mm(D) X 6 mm(P) X 4 mm(H)	15	ILOA003	CERTAIN LOCATOR ABUTMENT 4.1 mm(D) X 3 mm(H)	11
ART1178	PATIENT ED REGENERATIVE BROCHURE	33	ILOA004	CERTAIN LOCATOR ABUTMENT 4.1 mm(D) X 4 mm(H)	11
ART885	PATIENT ED BROCHURE - BONE LOSS (25-PACK)	33	ILOA005	CERTAIN LOCATOR ABUTMENT 4.1 mm(D) X 5 mm(H)	11
ART888	PATIENT ED BROCHURE - IMPLANT RETAINED, FIXED (25-PACK)	33	ILOA006	CERTAIN LOCATOR ABUTMENT 4.1 mm(D) X 6 mm(H)	11
ART889	PATIENT ED BROCHURE - IMPLANT RETAINED, LOCATOR® (25-PACK)	33	ILPAC3217	CERTAIN LOW PROFILE 17° PRE-ANGLED ABUTMENT 3.4 mm(D) X 2 mm(H)	9
ART890	PATIENT ED BROCHURE - IMPLANT RETAINED, BAR (25-PACK)	33	ILPAC3330	CERTAIN LOW PROFILE 30° PRE-ANGLED ABUTMENT 3.4 mm(D) X 3 mm(H)	9
ART929	PATIENT ED BROCHURE - IMMEDIATE TEETH REPLACEMENT, DIEM® (25-PACK)	33	ILPAC3417	CERTAIN LOW PROFILE 17° PRE-ANGLED ABUTMENT 3.4 mm(D) X 4 mm(H)	9
ART930	PATIENT ED CARE CARD (25-PACK)	33	ILPAC3530	CERTAIN LOW PROFILE 30° PRE-ANGLED ABUTMENT 3.4 mm(D) X 5 mm(H)	9
ART943	PATIENT ED FLIP CHART	33	ILPAC4217	CERTAIN LOW PROFILE 17° PRE-ANGLED ABUTMENT 4.1 mm(D) X 2 mm(H)	9
ART951EUR	PATIENT ED DVD - MULTILINGUAL VERSION	33	ILPAC4330	CERTAIN LOW PROFILE 30° PRE-ANGLED ABUTMENT 4.1 mm(D) X 3 mm(H)	9
ART957S	PATIENT ED BROCHURE - BONE LOSS (SPANISH VERSION)	33	ILPAC4417	CERTAIN LOW PROFILE 17° PRE-ANGLED ABUTMENT 4.1 mm(D) X 4 mm(H)	9
C9980	LOW TORQUE INDICATING RATCHET WRENCH STANDARD ISO 1797 ADAPTER	30	ILPAC4530	CERTAIN LOW PROFILE 30° PRE-ANGLED ABUTMENT 4.1 mm(D) X 5 mm(H)	9
C9981	LOW TORQUE INDICATING RATCHET WRENCH SHORT ISO 1797 ADAPTER	30	ILPAC5217	CERTAIN LOW PROFILE 17° PRE-ANGLED ABUTMENT 5 mm(D) X 2 mm(H)	9
CATDH	CONTRA ANGLE TORQUE DRIVER HANDLE	27	ILPAC5330	CERTAIN LOW PROFILE 30° PRE-ANGLED ABUTMENT 5 mm(D) X 3 mm(H)	9
DIEM2-4UNIT	PATIENT ED MODEL - DIEM 2 (4-UNIT)	33	ILPAC5417	CERTAIN LOW PROFILE 17° PRE-ANGLED ABUTMENT 5 mm(D) X 4 mm(H)	9
DIEM2-6UNIT	PATIENT ED MODEL - DIEM 2 (6-UNIT)	33	ILPAC5530	CERTAIN LOW PROFILE 30° PRE-ANGLED ABUTMENT 5 mm(D) X 5 mm(H)	9
EXHEXMOD	PATIENT ED MODEL - OSSEOTITE® IMPLANT	33	ILPC341	CERTAIN LOW PROFILE TWO-PIECE ABUTMENT 3.4 mm(D) X 1 mm(H)	9
GSH20	GOLD-TITE® (HEXED) RETAINING SCREW 2 mm	25,26	ILPC341U	CERTAIN LOW PROFILE ONE-PIECE ABUTMENT 3.4 mm(D) X 1 mm(H)	9
GSH30	GOLD-TITE (HEXED) RETAINING SCREW 3 mm	25,26	ILPC342	CERTAIN LOW PROFILE TWO-PIECE ABUTMENT 3.4 mm(D) X 2 mm(H)	9
GSH70	GOLD-TITE (HEXED) RETAINING SCREW 7 mm	25,26	ILPC342U	CERTAIN LOW PROFILE ONE-PIECE ABUTMENT 3.4 mm(D) X 2 mm(H)	9
GUCA1C	UCLA GOLD (HEXED) ABUTMENT CYLINDER 4.1 mm(D)	17,19	ILPC343	CERTAIN LOW PROFILE TWO-PIECE ABUTMENT 3.4 mm(D) X 3 mm(H)	9
GUCA2C	UCLA GOLD (NON-HEXED) ABUTMENT CYLINDER 4.1 mm(D)	19	ILPC343U	CERTAIN LOW PROFILE ONE-PIECE ABUTMENT 3.4 mm(D) X 3 mm(H)	9
HTD-C	CONTRA ANGLE TORQUE DRIVER BODY	27	ILPC344	CERTAIN LOW PROFILE TWO-PIECE ABUTMENT 3.4 mm(D) X 4 mm(H)	9
IAPP452G	CERTAIN® GINGIHUE POST 4.1 mm(D) X 5 mm(P) X 2 mm(H)	3	ILPC344U	CERTAIN LOW PROFILE ONE-PIECE ABUTMENT 3.4 mm(D) X 4 mm(H)	9
IAPP454G	CERTAIN GINGIHUE POST 4.1 mm(D) X 5 mm(P) X 4 mm(H)	3	ILPC441	CERTAIN LOW PROFILE TWO-PIECE ABUTMENT 4.1 mm(D) X 1 mm(H)	9
IAPP462G	CERTAIN GINGIHUE POST 4.1 mm(D) X 6 mm(P) X 2 mm(H)	3	ILPC441U	CERTAIN LOW PROFILE ONE-PIECE ABUTMENT 4.1 mm(D) X 1 mm(H)	9
IAPP464G	CERTAIN GINGIHUE POST 4.1 mm(D) X 6 mm(P) X 4 mm(H)	3	ILPC442	CERTAIN LOW PROFILE TWO-PIECE ABUTMENT 4.1 mm(D) X 2 mm(H)	9
ICD00	NARROW IMPRESSION COPING DRIVER	27,29	ILPC442U	CERTAIN LOW PROFILE ONE-PIECE ABUTMENT 4.1 mm(D) X 2 mm(H)	9
IGUCA1C	CERTAIN UCLA GOLD (HEXED) ABUTMENT CYLINDER 4.1 mm(D)	5,7	ILPC443	CERTAIN LOW PROFILE TWO-PIECE ABUTMENT 4.1 mm(D) X 3 mm(H)	9
IGUCA2C	CERTAIN UCLA GOLD (NON-HEXED) ABUTMENT CYLINDER 4.1 mm(D)	7	ILPC443U	CERTAIN LOW PROFILE ONE-PIECE ABUTMENT 4.1 mm(D) X 3 mm(H)	9
IGUCA2T	CERTAIN UCLA GOLD (NON-HEXED) ABUTMENT CYLINDER 4.1 mm(D) W/LARGE DIAMETER TITANIUM SCREW	7	ILPC444	CERTAIN LOW PROFILE TWO-PIECE ABUTMENT 4.1 mm(D) X 4 mm(H)	9
IIC12	IMPLANT EP® PICK-UP COPING 4.1 mm(D) X 5 mm(P)	16,18,20	ILPC444U	CERTAIN LOW PROFILE ONE-PIECE ABUTMENT 4.1 mm(D) X 4 mm(H)	9
IIC45	IMPLANT EP TWIST LOCK™ COPING 4.1 mm(D) X 5 mm(P)	16,18,20	ILPC541	CERTAIN LOW PROFILE TWO-PIECE ABUTMENT 5 mm(D) X 1 mm(H)	9
IIC46	IMPLANT EP TWIST LOCK COPING 4.1 mm(D) X 6 mm(P)	16,18,20	ILPC541U	CERTAIN LOW PROFILE ONE-PIECE ABUTMENT 5 mm(D) X 1 mm(H)	9
IIC60	IMPLANT EP PICK-UP COPING 4.1 mm(D) X 6 mm(P)	16,18,20	ILPC542	CERTAIN LOW PROFILE TWO-PIECE ABUTMENT 5 mm(D) X 2 mm(H)	9
IIC12	CERTAIN EP PICK-UP COPING 4.1 mm(D) X 5 mm(P)	4,6,8	ILPC542U	CERTAIN LOW PROFILE ONE-PIECE ABUTMENT 5 mm(D) X 2 mm(H)	9
IIC41	CERTAIN STRAIGHT PICK-UP COPING 4.1 mm(D) X 4.1 mm(P)	4,6,8	ILPC543	CERTAIN LOW PROFILE TWO-PIECE ABUTMENT 5 mm(D) X 3 mm(H)	9
IIC42	CERTAIN EP PICK-UP COPING 4.1 mm(D) X 4.1 mm(P)	4,6,8	ILPC543U	CERTAIN LOW PROFILE ONE-PIECE ABUTMENT 5 mm(D) X 3 mm(H)	9
			ILPC544	CERTAIN LOW PROFILE TWO-PIECE ABUTMENT 5 mm(D) X 4 mm(H)	9

Cat.#	Description	Page#	Cat.#	Description	Page#
ILPC544U	CERTAIN LOW PROFILE ONE-PIECE ABUTMENT 5 mm(D) X 4 mm(H)	9	IPAP562G	CERTAIN GINGIHUE 15° PRE-ANGLED POST 5 mm(D) X 6 mm(P) X 2 mm(H)	3
ILRGHG	CERTAIN GOLD-TITE (HEXED) LARGE SCREW	8,25	IPAP564G	CERTAIN GINGIHUE 15° PRE-ANGLED POST 5 mm(D) X 6 mm(P) X 4 mm(H)	3
ILRGHT	CERTAIN TITANIUM (HEXED) LARGE SCREW	8,25	IPAP662G	CERTAIN GINGIHUE 15° PRE-ANGLED POST 6 mm(D) X 6 mm(P) X 2 mm(H)	3
ILT150	CERTAIN LAPPING TOOL FOR CASTABLE UCLA ABUTMENT W/GUIDE 4 mm(D)	30	IPAP664G	CERTAIN GINGIHUE 15° PRE-ANGLED POST 6 mm(D) X 6 mm(P) X 4 mm(H)	3
ILTAH57	CERTAIN LABORATORY ABUTMENT HOLDER FOR ALL CERTAIN PLATFORMS	4,6,8,30	IPAPF454	CERTAIN PERFORMANCE 15° PRE-ANGLED POST 4.1 mm(D) X 5 mm(P) X 4 mm(H)	1
IMAP32G	CERTAIN GINGIHUE POST 3.4 mm (D) X 3.8 mm(P) X 2 mm(H)	3	IPAPF456	CERTAIN PERFORMANCE 15° PRE-ANGLED POST 4.1 mm(D) X 5 mm(P) X 6 mm(H)	1
IMAP34G	CERTAIN GINGIHUE POST 3.4 mm(D) X 3.8 mm(P) X 4 mm(H)	3	IPAPF554	CERTAIN PERFORMANCE 15° PRE-ANGLED POST 5 mm(D) X 5 mm(P) X 4 mm(H)	1
IMIC30	CERTAIN® EP® PICK-UP COPING 3.4 mm(D) X 3.4 mm(P)	4,6,8	IPAPF556	CERTAIN PERFORMANCE 15° PRE-ANGLED POST 5 mm(D) X 5 mm(P) X 6 mm(H)	1
IMIC33	CERTAIN EP PICK-UP COPING 3.4 mm(D) X 3.8 mm(P)	4,6,8	IPAPF664	CERTAIN® PERFORMANCE® 15° PRE-ANGLED POST 6 mm(D) X 6 mm(P) X 4 mm(H)	1
IMIC35	CERTAIN EP PICK-UP COPING 3.4 mm(D) X 5 mm(P)	4,6,8	IPAPF666	CERTAIN PERFORMANCE 15° PRE-ANGLED POST 6 mm(D) X 6 mm(P) X 6 mm(H)	1
IMIT33	CERTAIN EP TWIST LOCK COPING 3.4 mm(D) X 3.8 mm(P)	4,6,8	IPFP454	CERTAIN PERFORMANCE POST 4.1 mm(D) X 5 mm(P) X 4 mm(H)	1
IMIT35	CERTAIN EP TWIST LOCK COPING 3.4 mm(D) X 5 mm(P)	4,6,8	IPFP456	CERTAIN PERFORMANCE POST 4.1 mm(D) X 5 mm(P) X 6 mm(H)	1
IMLOA001	CERTAIN LOCATOR® ABUTMENT 3.4 mm(D) X 1 mm(H)	11	IPFTC41	CERTAIN PERFORMANCE (HEXED) TEMPORARY CYLINDER 4.1 mm(D)	2
IMLOA002	CERTAIN LOCATOR ABUTMENT 3.4 mm(D) X 2 mm(H)	11	IPFTC42	CERTAIN PERFORMANCE (NON-HEXED) TEMPORARY CYLINDER 4.1 mm(D)	2
IMLOA003	CERTAIN LOCATOR ABUTMENT 3.4 mm(D) X 3 mm(H)	11	IPPIA3	CERTAIN UCLA POLISHING PROTECTOR 4.1 mm(D)	30
IMLOA004	CERTAIN LOCATOR ABUTMENT 3.4 mm(D) X 4 mm(H)	11	IPPIA5	CERTAIN UCLA POLISHING PROTECTOR 5 mm(D)	30
IMLOA005	CERTAIN LOCATOR ABUTMENT 3.4 mm(D) X 5 mm(H)	11	IPPIA6	CERTAIN UCLA POLISHING PROTECTOR 6 mm(D)	30
IMLOA006	CERTAIN LOCATOR ABUTMENT 3.4 mm(D) X 6 mm(H)	11	IPPM1	CERTAIN UCLA POLISHING PROTECTOR 3.4 mm(D)	30
IMLT150	CERTAIN LAPPING TOOL FOR CASTABLE UCLA ABUTMENT W/GUIDE 3.4 mm(D)	31	IQSA01	CERTAIN QUICKSEAT® ACTIVATOR TOOL	7,30
IMMILA	CERTAIN IMPLANT LAB ANALOG 3.4 mm(D)	4,6,8	IRH600	CERTAIN REAMER WITH HANDLE	30
IMMTCS1	CERTAIN TEMPORARY HEXED CYLINDER 3.4 mm(D) W/HEXED SCREW	2	ISRT05	CERTAIN GUIDE HANDLE	30
IMPAP32G	CERTAIN GINGIHUE® 15° PRE-ANGLED POST 3.4 mm(D) X 3.8 mm(P) X 2 mm(H)	3	ISRT06	CERTAIN SCREW REMOVAL TOOL	30
IMPAP34G	CERTAIN GINGIHUE 15° PRE-ANGLED POST 3.4 mm(D) X 3.8 mm(P) X 4 mm(H)	3	ISRT07	MANUAL REVERSING TOOL	30
IMPAPF34	CERTAIN PERFORMANCE® 15° PRE-ANGLED POST 3.4 mm(D) X 3.8 mm(P) X 4 mm(H)	1	ITCS41	TEMPORARY (HEXED) IMPLANT CYLINDER 4.1 mm(D)	14
IMPAPF36	CERTAIN PERFORMANCE 15° PRE-ANGLED POST 3.4 mm(D) X 3.8 mm(P) X 6 mm(H)	1	ITCS42	TEMPORARY (NON-HEXED) IMPLANT CYLINDER 4.1 mm(D)	14
IMPFP34	CERTAIN PERFORMANCE POST 3.4 mm(D) X 3.8 mm(P) X 4 mm(H)	1	IUNAB2C	CERTAIN UCLA CASTABLE (NON-HEXED) ABUTMENT CYLINDER 4.1 mm(D)	7
IMPFP36	CERTAIN PERFORMANCE POST 3.4 mm(D) X 3.8 mm(P) X 6 mm(H)	1	IUNAB2T	CERTAIN UCLA CASTABLE (NON-HEXED) ABUTMENT CYLINDER 4.1 mm(D) W/LARGE DIAMETER TITANIUM SCREW	7
IMPFTC32	CERTAIN PERFORMANCE (NON-HEXED) TEMPORARY CYLINDER 3.4 mm(D)	2	IUNIHG	CERTAIN GOLD-TITE® (HEXED) SCREW	4,6,8,25
IMPFTC34	CERTAIN PERFORMANCE (HEXED) TEMPORARY CYLINDER 3.4 mm(D)	2	IUNIHT	CERTAIN TITANIUM (HEXED) SCREW	25
IMUCC2C	CERTAIN UCLA CASTABLE (NON-HEXED) ABUTMENT CYLINDER 3.4 mm(D)	7	IUNITS	CERTAIN TRY-IN (HEXED) SCREW (5-PACK) 3.4, 4.1, 5, 6 mm(D)	4,6,8,25
IMUCC2T	CERTAIN UCLA CASTABLE (NON-HEXED) ABUTMENT CYLINDER 3.4 mm(D) W/LARGE DIAMETER TITANIUM SCREW	7	IWGA51C	CERTAIN UCLA GOLD (HEXED) ABUTMENT CYLINDER 5 mm(D)	5,7
IMUCG1C	CERTAIN UCLA GOLD (HEXED) ABUTMENT CYLINDER 3.4 mm(D)	5,7	IWGA52C	CERTAIN UCLA GOLD (NON-HEXED) ABUTMENT CYLINDER 5 mm(D)	7
IMUCG2C	CERTAIN UCLA GOLD (NON-HEXED) ABUTMENT CYLINDER 3.4 mm(D)	7	IWGA52T	CERTAIN UCLA GOLD (NON-HEXED) ABUTMENT CYLINDER 5 mm(D) W/LARGE DIAMETER TITANIUM SCREW	7
IMUCG2T	CERTAIN UCLA GOLD (NON-HEXED) ABUTMENT CYLINDER 3.4 mm(D) W/LARGE DIAMETER TITANIUM SCREW	7	IWGA61C	CERTAIN UCLA GOLD (HEXED) ABUTMENT CYLINDER 6 mm(D)	5,7
IPAP452G	CERTAIN GINGIHUE 15° PRE-ANGLED POST 4.1 mm(D) X 5 mm(P) X 2 mm(H)	3	IWGA62C	CERTAIN UCLA GOLD (NON-HEXED) ABUTMENT CYLINDER 6 mm(D)	7
IPAP454G	CERTAIN GINGIHUE 15° PRE-ANGLED POST 4.1 mm(D) X 5 mm(P) X 4 mm(H)	3	IWGA62T	CERTAIN UCLA GOLD (NON-HEXED) ABUTMENT CYLINDER 6 mm(D) W/LARGE DIAMETER TITANIUM SCREW	7
IPAP462G	CERTAIN GINGIHUE 15° PRE-ANGLED POST 4.1 mm(D) X 6 mm(P) X 2 mm(H)	3	IWIP50	CERTAIN STRAIGHT PICK-UP COPING 5 mm(D) X 5 mm(P)	4,6,8
IPAP464G	CERTAIN GINGIHUE 15° PRE-ANGLED POST 4.1 mm(D) X 6 mm(P) X 4 mm(H) 4.1 mm(D) X 7.5 mm(P) X 4 mm(H)	3	IWIP55	CERTAIN EP® PICK-UP COPING 5 mm(D) X 5.6 mm(P)	4,6,8
IPAP552G	CERTAIN GINGIHUE 15° PRE-ANGLED POST 5 mm(D) X 5 mm(P) X 2 mm(H)	3	IWIP56	CERTAIN EP PICK-UP COPING 5 mm(D) X 6 mm(P)	4,6,8
IPAP554G	CERTAIN GINGIHUE 15° PRE-ANGLED POST 5 mm(D) X 5 mm(P) X 4 mm(H)	3	IWIP60	CERTAIN STRAIGHT PICK-UP COPING 6 mm(D) X 6 mm(P)	4,6,8
			IWIP66	CERTAIN EP PICK-UP COPING 6 mm(D) X 6.6 mm(P)	4,6,8
			IWIT50	CERTAIN STRAIGHT TRANSFER COPING 5 mm(D) X 5 mm(P)	4,6,8
			IWIT55	CERTAIN EP TWIST LOCK™ COPING 5 mm(D) X 5.6 mm(P)	4,6,8
			IWIT56	CERTAIN EP TWIST LOCK COPING 5 mm(D) X 6 mm(P)	4,6,8

Index

Cat.#	Description	Page#	Cat.#	Description	Page#
IWIT60	CERTAIN STRAIGHT TRANSFER COPING 6 mm(D) X 6 mm(P)	4,6,8	LPAC4217	LOW PROFILE 17° PRE-ANGLED ABUTMENT 4.1 mm(D) X 2 mm(H)	21
IWIT66	CERTAIN EP TWIST LOCK COPING 6 mm(D) X 6.6 mm(P)	4,6,8	LPAC4330	LOW PROFILE 30° PRE-ANGLED ABUTMENT 4.1 mm(D) X 3 mm(H)	21
IWPC52C	CERTAIN UCLA CASTABLE (NON-HEXED) ABUTMENT CYLINDER 5 mm(D)	7	LPAC4417	LOW PROFILE 17° PRE-ANGLED ABUTMENT 4.1 mm(D) X 4 mm(H)	21
IWPC52T	CERTAIN UCLA CASTABLE (NON-HEXED) ABUTMENT CYLINDER 5 mm(D) W/LARGE DIAMETER TITANIUM SCREW	7	LPAC4530	LOW PROFILE 30° PRE-ANGLED ABUTMENT 4.1 mm(D) X 5 mm(H)	21
IWPC62C	CERTAIN UCLA CASTABLE (NON-HEXED) ABUTMENT CYLINDER 6 mm(D)	7	LPAC5217	LOW PROFILE 17° PRE-ANGLED ABUTMENT 5 mm(D) X 2 mm(H)	21
IWPC62T	CERTAIN UCLA CASTABLE (NON-HEXED) ABUTMENT CYLINDER 6 mm(D) W/LARGE DIAMETER TITANIUM SCREW	7	LPAC5330	LOW PROFILE 30° PRE-ANGLED ABUTMENT 5 mm(D) X 3 mm(H)	21
IWPPF554	CERTAIN PREFORMANCE POST 5 mm(D) X 5 mm(P) X 4 mm(H)	1	LPAC5417	LOW PROFILE 17° PRE-ANGLED ABUTMENT 5 mm(D) X 4 mm(H)	21
IWPPF556	CERTAIN PREFORMANCE POST 5 mm(D) X 5 mm(P) X 6 mm(H)	1	LPAC5530	LOW PROFILE 30° PRE-ANGLED ABUTMENT 5 mm(D) X 5 mm(H)	21
IWPPF664	CERTAIN PREFORMANCE POST 6 mm(D) X 6 mm(P) X 4 mm(H)	1	LPC341	LOW PROFILE TWO-PIECE ABUTMENT 3.4 mm(D) X 1 mm(H)	21
IWPPF666	CERTAIN PREFORMANCE POST 6 mm(D) X 6 mm(P) X 6 mm(H)	1	LPC341U	LOW PROFILE ONE-PIECE ABUTMENT 3.4 mm(D) X 1 mm(H)	21
IWPFTC51	CERTAIN PREFORMANCE (HEXED) TEMPORARY CYLINDER 5 mm(D)	2	LPC342	LOW PROFILE TWO-PIECE ABUTMENT 3.4 mm(D) X 2 mm(H)	21
IWPFTC52	CERTAIN PREFORMANCE (NON-HEXED) TEMPORARY CYLINDER 5 mm(D)	2	LPC342U	LOW PROFILE ONE-PIECE ABUTMENT 3.4 mm(D) X 2 mm(H)	21
IWPFTC61	CERTAIN PREFORMANCE (HEXED) TEMPORARY CYLINDER 6 mm(D)	2	LPC343	LOW PROFILE TWO-PIECE ABUTMENT 3.4 mm(D) X 3 mm(H)	21
IWPFTC62	CERTAIN PREFORMANCE (NON-HEXED) TEMPORARY CYLINDER 6 mm(D)	2	LPC343U	LOW PROFILE ONE-PIECE ABUTMENT 3.4 mm(D) X 3 mm(H)	21
IWPP552G	CERTAIN GINGIHUE® POST 5 mm(D) X 5 mm(P) X 2 mm(H)	3	LPC344	LOW PROFILE TWO-PIECE ABUTMENT 3.4 mm(D) X 4 mm(H)	21
IWPP554G	CERTAIN GINGIHUE POST 5 mm(D) X 5 mm(P) X 4 mm(H)	3	LPC344U	LOW PROFILE ONE-PIECE ABUTMENT 3.4 mm(D) X 4 mm(H)	21
IWPP562G	CERTAIN GINGIHUE POST 5 mm(D) X 6 mm(P) X 2 mm(H)	3	LPC441	LOW PROFILE TWO-PIECE ABUTMENT 4.1 mm(D) X 1 mm(H)	21
IWPP564G	CERTAIN GINGIHUE POST 5 mm(D) X 6 mm(P) X 4 mm(H)	3	LPC441U	LOW PROFILE ONE-PIECE ABUTMENT 4.1 mm(D) X 1 mm(H)	21
IWPP662G	CERTAIN GINGIHUE POST 6 mm(D) X 6 mm(P) X 2 mm(H)	3	LPC442	LOW PROFILE TWO-PIECE ABUTMENT 4.1 mm(D) X 2 mm(H)	21
IWPP664G	CERTAIN GINGIHUE POST 6 mm(D) X 6 mm(P) X 4 mm(H)	3	LPC442U	LOW PROFILE ONE-PIECE ABUTMENT 4.1 mm(D) X 2 mm(H)	21
IWSU30	CERTAIN WAXING SCREW W/GUIDE PIN	25,30	LPC443	LOW PROFILE TWO-PIECE ABUTMENT 4.1 mm(D) X 3 mm(H)	21
IWTCS51	CERTAIN TEMPORARY (HEXED) CYLINDER 5 mm(D) W/HEXED SCREW	2	LPC443U	LOW PROFILE ONE-PIECE ABUTMENT 4.1 mm(D) X 3 mm(H)	21
IWTCS52	CERTAIN TEMPORARY (NON-HEXED) CYLINDER 5 mm(D) W/LARGE DIA. SCREW	2	LPC444	LOW PROFILE TWO-PIECE ABUTMENT 4.1 mm(D) X 4 mm(H)	21
IWTCS61	CERTAIN TEMPORARY (HEXED) CYLINDER 6 mm(D) W/HEXED SCREW	2	LPC444U	LOW PROFILE ONE-PIECE ABUTMENT 4.1 mm(D) X 4 mm(H)	21
IWTCS62	CERTAIN TEMPORARY (NON-HEXED) CYLINDER 6 mm(D) W/LARGE DIA. SCREW	2	LPC541	LOW PROFILE TWO-PIECE ABUTMENT 5 mm(D) X 1 mm(H)	21
LAELM	LOCATOR® EXTRA LIGHT RETENTION MALES (4-PACK)	12,24	LPC542	LOW PROFILE TWO-PIECE ABUTMENT 5 mm(D) X 2 mm(H)	21
LAERM	LOCATOR EXTENDED RANGE MALES - GREEN, 4LBS/1800G (4-PACK)	12,24	LPC542U	LOW PROFILE ONE-PIECE ABUTMENT 5 mm(D) X 2 mm(H)	21
LAIC1	LOCATOR ABUTMENT IMPRESSION COPING	12,24	LPC543	LOW PROFILE TWO-PIECE ABUTMENT 5 mm(D) X 3 mm(H)	21
LALA1	LOCATOR ABUTMENT LAB ANALOG	12,24	LPC543U	LOW PROFILE ONE-PIECE ABUTMENT 5 mm(D) X 3 mm(H)	21
LARMS	LOCATOR REPLACEMENT MALES - WHITE, 5LBS/2250G (4-PACK)	12,24	LPC544	LOW PROFILE TWO-PIECE ABUTMENT 5 mm(D) X 4 mm(H)	21
LCTDR1	LOCATOR CORE TOOL/ABUTMENT DRIVER	12,24	LPC544U	LOW PROFILE ONE-PIECE ABUTMENT 5 mm(D) X 4 mm(H)	21
LELARM	EXTRA LIGHT RETENTION EXTENDED RANGE MALES - RED, 1.5LBS/680G	12,24	LPCAMI	LOW PROFILE ABUTMENT LAPPING TOOL	10,22
LLRMS	LOCATOR LIGHT RETENTION REPLACEMENT MALES - PINK, 3LBS/1350G (4-PACK)	12,24	LPCC1	LOW PROFILE ABUTMENT (HEXED) CASTABLE CYLINDER	10,22
LOA001	LOCATOR ABUTMENT 4.1 mm(D) X 1 mm(H)	23	LPCC2	LOW PROFILE ABUTMENT (NON-HEXED) CASTABLE CYLINDER	10,22
LOA002	LOCATOR ABUTMENT 4.1 mm(D) X 2 mm(H)	23	LPCDE	LOW PROFILE ABUTMENT DISTAL EXTENSION	10,22
LOA003	LOCATOR ABUTMENT 4.1 mm(D) X 3 mm(H)	23	LPCGC1	LOW PROFILE ABUTMENT (HEXED) GOLD CYLINDER	10,22
LOA004	LOCATOR ABUTMENT 4.1 mm(D) X 4 mm(H)	23	LPCGC2	LOW PROFILE ABUTMENT (NON-HEXED) GOLD CYLINDER	10,22
LOA005	LOCATOR ABUTMENT 4.1 mm(D) X 5 mm(H)	23	LPCGSH	LOW PROFILE ABUTMENT GOLD-TITE® RETAINING SCREW	10,22,25,26
LOA006	LOCATOR ABUTMENT 4.1 mm(D) X 6 mm(H)	23	LPCTSH	LOW PROFILE ABUTMENT TITANIUM RETAINING SCREW	10,22,25,26
LOAB	LOCATOR BAR ATTACHMENT	12,24	LPCHC	LOW PROFILE ABUTMENT HEALING CAP	10,22
LOAD	LOCATOR BAR ATTACHMENT DRILL	12,24	LPCLA	LOW PROFILE ABUTMENT LAB ANALOG	10,22
LOADT4	LOCATOR ABUTMENT DRIVER TIP 24 mm(L)	12,24,29	LPCPIC1	LOW PROFILE ABUTMENT (HEXED) PICK-UP IMPRESSION COPING	10,22
LOADT9	LOCATOR ABUTMENT DRIVER TIP 30 mm(L)	12,24,29	LPCPIC2	LOW PROFILE ABUTMENT (NON-HEXED) PICK-UP IMPRESSION COPING	10,22
LOAH	LOCATOR REPLACEMENT HOUSING	12,24	LPCPP	LOW PROFILE ABUTMENT POLISHING PROTECTOR	10,22
LOAT	LOCATOR BAR ATTACHMENT TAP	12,24	LPCPTC1	LOW PROFILE ABUTMENT (HEXED) PREFORMANCE® TEMPORARY CYLINDER	10,22
LORHK	LOCATOR REPLACEMENT HOUSING KIT	12,24	LPCPTC2	LOW PROFILE ABUTMENT (NON-HEXED) PREFORMANCE TEMPORARY CYLINDER	10,22
LPAC3217	LOW PROFILE 17° PRE-ANGLED ABUTMENT 3.4 mm(D) X 2 mm(H)	21	LPCQB	LOW PROFILE ABUTMENT QUICKBRIDGE®	10,22
LPAC3330	LOW PROFILE 30° PRE-ANGLED ABUTMENT 3.4 mm(D) X 3 mm(H)	21	LPCQBCAP	LOW PROFILE ABUTMENT QUICKBRIDGE REPLACEMENT CAP	10,22
LPAC3417	LOW PROFILE 17° PRE-ANGLED ABUTMENT 3.4 mm(D) X 4 mm(H)	21	LPCRIC	LOW PROFILE ABUTMENT QUICKBRIDGE IMPRESSION COPING	10,22
LPAC3530	LOW PROFILE 30° PRE-ANGLED ABUTMENT 3.4 mm(D) X 5 mm(H)	21	LPCTC1	LOW PROFILE ABUTMENT (HEXED) TEMPORARY CYLINDER	10,22
			LPCTC2	LOW PROFILE ABUTMENT (NON-HEXED) TEMPORARY CYLINDER	10,22
			LPCTIC1	LOW PROFILE ABUTMENT (HEXED) TRANSFER IMPRESSION COPING	10,22

Cat.#	Description	Page#	Cat.#	Description	Page#
LPCTIC2	LOW PROFILE ABUTMENT (NON-HEXED) TRANSFER IMPRESSION COPING	10,22	PAP564G	GINGIHUE 15° PRE-ANGLED POST 5 mm(D) X 6 mm(P) X 4 mm(H)	15
LPCTSH	LOW PROFILE ABUTMENT TITANIUM RETAINING SCREW	10,22	PAP662G	GINGIHUE 15° PRE-ANGLED POST 6 mm(D) X 6 mm(P) X 2 mm(H)	15
LPCTUCA	LOW PROFILE ABUTMENT LASER WELDED CYLINDER	10,22	PAP664G	GINGIHUE 15° PRE-ANGLED POST 6 mm(D) X 6 mm(P) X 4 mm(H)	15
LPCWS	LOW PROFILE ABUTMENT WAXING SCREW	10,22	PAPF454	PREFORMANCE® 13° PRE-ANGLED POST 4.1 mm(D) X 5 mm(P) X 4 mm(H)	13
LT150	LAPPING TOOL FOR CASTABLE UCLA ABUTMENT W/GUIDE 4 mm(D)	30	PAPF456	PREFORMANCE 13° PRE-ANGLED POST 4.1 mm(D) X 5 mm(P) X 6 mm(H)	13
LTAH5	LABORATORY ABUTMENT HOLDER FOR 3.4 mm(D)	16,18,20,30	PAPF554	PREFORMANCE 13° PRE-ANGLED POST 5 mm(D) X 5 mm(P) X 4 mm(H)	13
LTAH7	LABORATORY HOLDER FOR 4.1, 5, 6 mm(D)	16,18,20,30	PAPF556	PREFORMANCE 13° PRE-ANGLED POST 5 mm(D) X 5 mm(P) X 6 mm(H)	13
L-TIRW	LOW TORQUE INDICATING RATCHET WRENCH	27,30	PAPF664	PREFORMANCE 13° PRE-ANGLED POST 6 mm(D) X 6 mm(P) X 4 mm(H)	13
MAP32G	GINGIHUE® POST 3.4 mm(D) X 3.8 mm(P) X 2 mm(H)	15	PAPF666	PREFORMANCE 13° PRE-ANGLED POST 6 mm(D) X 6 mm(P) X 6 mm(H)	13
MAP34G	GINGIHUE POST 3.4 mm(D) X 3.8 mm(P) X 4 mm(H)	15	PEMODLOC	PATIENT ED MODEL - LOCATOR®	33
MIC33	PICK-UP IMPRESSION COPING 3.4 mm(D)	16,18,20	PEMODSM	PATIENT ED MODEL - SINGLE/MULTIPLE TEETH	33
MIDTH	MANUAL INTERCHANGEABLE DRIVER TIP HANDLE	29	PFP454	PREFORMANCE POST 4.1 mm(D) X 5 mm(P) X 4 mm(H)	13
MLOA001	LOCATOR® ABUTMENT 3.4 mm(D) X 1 mm(H)	23	PFP456	PREFORMANCE POST 4.1 mm(D) X 5 mm(P) X 6 mm(H)	13
MLOA002	LOCATOR ABUTMENT 3.4 mm(D) X 2 mm(H)	23	PFTC41	PREFORMANCE (HEXED) TEMPORARY CYLINDER 4.1 mm(D)	14
MLOA003	LOCATOR ABUTMENT 3.4 mm(D) X 3 mm(H)	23	PFTC42	PREFORMANCE (NON-HEXED) TEMPORARY CYLINDER 4.1 mm(D)	14
MLOA004	LOCATOR ABUTMENT 3.4 mm(D) X 4 mm(H)	23	PHD00N	NARROW POSTERIOR SMALL HEX DRIVER 17 mm(L)	29
MLOA005	LOCATOR ABUTMENT 3.4 mm(D) X 5 mm(H)	23	PHD01N	NARROW STANDARD SMALL HEX DRIVER 24 mm(L)	29
MLOA006	LOCATOR ABUTMENT 3.4 mm(D) X 6 mm(H)	23	PHD02N	NARROW POSTERIOR LARGE HEX DRIVER 17 mm(L)	25,26,27,29
MLT150	LAPPING TOOL FOR CASTABLE UCLA ABUTMENT W/GUIDE 3.4 mm(D)	30	PHD03N	NARROW STANDARD LARGE HEX DRIVER 24 mm(L)	25,26,27,29
MMILA	LAB ANALOG 3.4 mm(D)	16,18,20	PPIA3	POLISHING PROTECTOR UCLA ABUTMENT 4.1 mm(D)	30
MMTCS1	TEMPORARY (HEXED) CYLINDER 3.4 mm(D) W/RETENTION AND (HEXED) SCREW	14	PPMM1	POLISHING PROTECTOR UCLA 3.4 mm(D)	30
MPAP32G	GINGIHUE 15° PRE-ANGLED POST 3.4 mm(D) X 3.8 mm(P) X 2 mm(H)	15	PSD00	POSTERIOR SLOTTED SCREW DRIVER - SLOTTED 17 mm(L)	27
MPAP34G	GINGIHUE 15° PRE-ANGLED POST 3.4 mm(D) X 3.8 mm(P) X 4 mm(H)	15	PSD01	STANDARD SLOTTED SCREW DRIVER - SLOTTED 24 mm(L)	27
MPAPF34	PREFORMANCE 15° PRE-ANGLED POST 3.4 mm(D) X 3.8 mm(P) X 4 mm(H)	13	PSDT1	PROSTHETIC SYSTEM DRIVER TRAY	27
MPAPF36	PREFORMANCE 15° PRE-ANGLED POST 3.4 mm(D) X 3.8 mm(P) X 6 mm(H)	13	PSQD0N	NARROW POSTERIOR SQUARE DRIVER 17 mm(L)	26,27,29
MPFP34	PREFORMANCE POST 3.4 mm(D) X 3.8 mm(P) X 4 mm(H)	13	PSQD1N	NARROW STANDARD SLOTTED SQUARE DRIVER 24 mm(L)	26,27,29
MPFP36	PREFORMANCE POST 3.4 mm(D) X 3.8 mm(P) X 6 mm(H)	13	RASA3	RIGHT ANGLE ABUTMENT DRIVER TIP - STEEL	27,29
MPFTC32	PREFORMANCE (NON-HEXED) TEMPORARY CYLINDER 3.4 mm(D)	14	RASD1	RIGHT ANGLE SLOTTED DRIVER TIP 24 mm(L)	27
MPFTC34	PREFORMANCE (HEXED) TEMPORARY CYLINDER 3.4 mm(D)	14	RASD6	RIGHT ANGLE SLOTTED DRIVER TIP 30 mm(L)	27
MSGIAM1	NAVIGATOR® CERTAIN® ANALOG MOUNT 3.4 mm(D) X 1(L)	31	RASH2N	NARROW RIGHT ANGLE SMALL (HEXED) DRIVER TIP 24 mm(L)	29
MSGIAM2	NAVIGATOR CERTAIN ANALOG MOUNT 3.4 mm(D) X 2(L)	31	RASH3N	NARROW RIGHT ANGLE LARGE (HEXED) DRIVER TIP 24 mm(L)	25,26,27,29
MSGIAM3	NAVIGATOR CERTAIN ANALOG MOUNT 3.4 mm(D) X 3(L)	31	RASH4	RIGHT ANGLE DRIVER TIP 50 IN. (HEXED) 24 mm(L)	12,24,29
MSGIAM4	NAVIGATOR CERTAIN ANALOG MOUNT 3.4 mm(D) X 4(L)	31	RASH7N	NARROW RIGHT ANGLE SMALL (HEXED) DRIVER TIP 30 mm(L)	29
MUCC1C	UCLA CASTABLE (HEXED) ABUTMENT CYLINDER 3.4 mm(D)	19	RASH8N	NARROW RIGHT ANGLE LARGE (HEXED) DRIVER TIP 30 mm(L)	25,26,27,29
MUCC2C	UCLA CASTABLE (NON-HEXED) ABUTMENT CYLINDER 3.4 mm(D)	19	RASH9	RIGHT ANGLE DRIVER TIP 50 IN. (HEXED) 30 mm(L)	12,24,29
MUCG1C	UCLA GOLD (HEXED) ABUTMENT CYLINDER 3.4 mm(D)	17,19	RASQ3N	NARROW RIGHT ANGLE SQUARE DRIVER TIP 24 mm(L)	26,27,29
MUCG2C	UCLA GOLD (NON-HEXED) ABUTMENT CYLINDER 3.4 mm(D)	19	RASQ8N	NARROW RIGHT ANGLE SQUARE DRIVER TIP 30 mm(L)	26,27,29
MUNITS	SQUARE TRY-IN SCREW (5-PACK), 3.4 mm(D)	16,18,20,26	RH600	REAMER/HANDLE FOR SGC34, UNAB1/2	30
NCATD0	CONTRA ANGLE TORQUE DRIVER KIT	27	RTI2035TR	RESTORATIVE TORQUE INDICATOR TRAY	27
NCATD0C	CONTRA ANGLE TORQUE DRIVER KIT CERTAIN	27	SGIAM41	NAVIGATOR® CERTAIN® ANALOG MOUNT 4.1 mm(D) X 1(L)	31
NPSDK0	PROSTHETIC DRIVER KIT	27	SGIAM42	NAVIGATOR CERTAIN ANALOG MOUNT 4.1 mm(D) X 2(L)	31
PAD00	POSTERIOR ABUTMENT DRIVER 17 mm(L)	27,29	SGIAM43	NAVIGATOR CERTAIN ANALOG MOUNT 4.1 mm(D) X 3(L)	31
PAD24	STANDARD ABUTMENT DRIVER 24 mm(L)	27,29	SGIAM44	NAVIGATOR CERTAIN ANALOG MOUNT 4.1 mm(D) X 4(L)	31
PAP452G	GINGIHUE® 15° PRE-ANGLED POST 4.1 mm(D) X 5 mm(P) X 2 mm(H)	15	SGIAM51	NAVIGATOR CERTAIN ANALOG MOUNT 5 mm(D) X 1(L)	31
PAP454G	GINGIHUE 15° PRE-ANGLED POST 4.1 mm(D) X 5 mm(P) X 4 mm(H)	15	SGIAM52	NAVIGATOR CERTAIN ANALOG MOUNT 5 mm(D) X 2(L)	31
PAP462G	GINGIHUE 15° PRE-ANGLED POST 4.1 mm(D) X 6 mm(P) X 2 mm(H)	15	SGIAM53	NAVIGATOR CERTAIN ANALOG MOUNT 5 mm(D) X 3(L)	31
PAP464G	GINGIHUE 15° PRE-ANGLED POST 4.1 mm(D) X 6 mm(P) X 4 mm(H)	15	SGIAM54	NAVIGATOR CERTAIN ANALOG MOUNT 5 mm(D) X 4(L)	31
PAP552G	GINGIHUE 15° PRE-ANGLED POST 5 mm(D) X 5 mm(P) X 2 mm(H)	15	SGLKIT	NAVIGATOR LABORATORY KIT	31
PAP554G	GINGIHUE 15° PRE-ANGLED POST 5 mm(D) X 5 mm(P) X 4 mm(H)	15	SGLTRAY	NAVIGATOR LABORATORY TRAY	31
PAP562G	GINGIHUE 15° PRE-ANGLED POST 5 mm(D) X 6 mm(P) X 2 mm(H)	15	SGUCA1C	UCLA GOLD STANDARD ZR (HEXED) ABUTMENT CYLINDER 4.1 mm(D)	17,19
			SRT02	HAND DRILL 14 mm(D) BRONZE DRILL FOR SRT10	30
			SRT03	HAND DRILL 1.18 mm(D) SILVER DRILL FOR SRT10	30
			SRT04	EXTRACTING DRILL - ISO LATCH FOR SRT10	30
			SRT05	GUIDE HANDLE FOR SRT10	30

Index

Cat.#	Description	Page#			
SWGA51C	UCLA GOLD STANDARD ZR (HEXED) ABUTMENT CYLINDER 5 mm(D)	17,19	WTCS61	TEMPORARY (HEXED) IMPLANT CYLINDER 6 mm(D)	14
SWGA61C	UCLA GOLD STANDARD ZR (HEXED) ABUTMENT CYLINDER 6 mm(D)	17,19	WTCS62	TEMPORARY (NON-HEXED) IMPLANT CYLINDER 6 mm(D)	14
TD-10	CONTRA ANGLE TORQUE CONTROLLER 10 NCM	27			
TD-20	CONTRA ANGLE TORQUE CONTROLLER 20 NCM	27			
TD-35	CONTRA ANGLE TORQUE CONTROLLER 32 NCM	27			
THRC4	TEMPORARY HEALING RETENTION CYLINDER 4.1 mm(D) X 4 mm(H)	23			
THRC6	TEMPORARY HEALING RETENTION CYLINDER 4.1 mm(D) X 6 mm(H)	23			
TMP80	TISSUE MEASURING POST	30			
UNA125	UCLA CASTABLE (HEXED) ABUTMENT CYLINDER (25-PACK)	17,19			
UNA225	UCLA CASTABLE (NON-HEXED) ABUTMENT CYLINDER (25-PACK)	19			
UNAB1C	UCLA CASTABLE (HEXED) ABUTMENT CYLINDER	17,19			
UNAB2C	UCLA CASTABLE (NON-HEXED) ABUTMENT CYLINDER	19			
UNIHG	GOLD-TITE® (HEXED) SCREW	16,18,20,26			
UNIHT	TITANIUM (HEXED) SCREW	16,18,20,26			
UNISG	GOLD-TITE SQUARE SCREW	16,18,20,26			
UNITS	SQUARE TRY-IN SCREW (5-PACK) 4.1, 5, 6 mm(D)	16,18,20,26			
WGA51C	UCLA GOLD (HEXED) ABUTMENT CYLINDER 5 mm(D)	17,19			
WGA52C	UCLA GOLD (NON-HEXED) ABUTMENT CYLINDER 5 mm(D)	19			
WGA61C	UCLA GOLD (HEXED) ABUTMENT CYLINDER 6 mm(D)	17,19			
WGA62C	UCLA GOLD (NON-HEXED) ABUTMENT CYLINDER 6 mm(D)	19			
WIP55	IMPLANT EP® PICK-UP COPING 5 mm(D) X 5 mm(P)	16,18,20			
WIP56	IMPLANT EP PICK-UP COPING 5 mm(D) X 6 mm(P)	16,18,20			
WIP66	IMPLANT EP PICK-UP COPING 6 mm(D) X 6 mm(P)	16,18,20			
WIT55	IMPLANT EP TWIST LOCK™ COPING 5 mm(D) X 5 mm(P)	16,18,20			
WIT56	IMPLANT EP TWIST LOCK COPING 5 mm(D) X 6 mm(P)	16,18,20			
WIT66	IMPLANT EP TWIST LOCK COPING 6 mm(D) X 6 mm(P)	16,18,20			
WLT10	LAPPING TOOL FOR WIDE DIAMETER IMPLANTS	30			
WPC51C	UCLA CASTABLE (HEXED) ABUTMENT CYLINDER 5 mm(D)	17,19			
WPC52C	UCLA CASTABLE (NON-HEXED) ABUTMENT CYLINDER 5 mm(D)	19			
WPC61C	UCLA CASTABLE (HEXED) ABUTMENT CYLINDER 6 mm(D)	17,19			
WPC62C	UCLA CASTABLE (NON-HEXED) ABUTMENT CYLINDER 6 mm(D)	19			
WPPF554	PREFORMANCE® POST 5 mm(D) X 5 mm(P) X 4 mm(H)	13			
WPPF556	PREFORMANCE POST 5 mm(D) X 5 mm(P) X 6 mm(H)	13			
WPPF664	PREFORMANCE POST 6 mm(D) X 6 mm(P) X 4 mm(H)	13			
WPPF666	PREFORMANCE POST 6 mm(D) X 6 mm(P) X 6 mm(H)	13			
WPFTC51	PREFORMANCE (HEXED) TEMPORARY CYLINDER 5 mm	14			
WPFTC52	PREFORMANCE (NON-HEXED) TEMPORARY CYLINDER 5 mm	14			
WPFTC61	PREFORMANCE (HEXED) TEMPORARY CYLINDER 6 mm	14			
WPFTC62	PREFORMANCE (NON-HEXED) TEMPORARY CYLINDER 6 mm	14			
WPP50	POLISHING PROTECTOR 5 mm(D)	31			
WPP552G	GINGIHUE® POST 5 mm(D) X 5 mm(P) X 2 mm(H)	15			
WPP554G	GINGIHUE POST 5 mm(D) X 5 mm(P) X 4 mm(H)	15			
WPP562G	GINGIHUE POST 5 mm(D) X 6 mm(P) X 2 mm(H)	15			
WPP564G	GINGIHUE POST 5 mm(D) X 6 mm(P) X 4 mm(H)	15			
WPP60	POLISHING PROTECTOR 6 mm(D)	30			
WPP662G	GINGIHUE POST 6 mm(D) X 6 mm(P) X 2 mm(H)	15			
WPP664G	GINGIHUE POST 6 mm(D) X 6 mm(P) X 4 mm(H)	15			
WSK10	ABUTMENT WAXING SCREW/GUIDE PIN-KNURLED 10 mm(L)	26			
WSK15	ABUTMENT WAXING SCREW/GUIDE PIN-KNURLED 15 mm(L)	26			
WSK30	ABUTMENT WAXING SCREW/GUIDE PIN-KNURLED 20 mm(L)	26			
WSU30	LABORATORY WAXING SCREW/GUIDE PIN-KNURLED	26			
WTCS51	TEMPORARY (HEXED) IMPLANT CYLINDER 5 mm(D)	14			
WTCS52	TEMPORARY (NON-HEXED) IMPLANT CYLINDER 5 mm(D)	14			

Ordering Form

Customer Account#

Date: ____/____/____

Representative taking order

Customer Name:

Specialty:

Bill To: _____

Phone #: () _____

Ship To:_____

Phone #: () _____

Would you like to receive information via email about
BIOMET **3i**'s Products, Services and Events?

☐ Yes ☐ No

Email Address: _____

Special Instructions:

Contact Name:

Customer P.O. #: _____

Method of Shipment: _____FedEx _____UPS

____Next Day by 8:30am

____Next Day by 10:30am

____Next Day - Afternoon

_____Two Day

_____Ground

_____ Other _____

Method of Payment:

☐ Credit Card: _____ Visa _____ MasterCard _____ AMEX
Card #: _____

Expiration Date: _____

☐ C.O.D.☐ Bill - Net 30[illegible]

Date:

Order Number:

Order Total:

Customer Service Representative:

BIOMET **3i**, LLC
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: 800-342-5454 • Fax: 561-776-1272
Website: www.biomet3i.com

Thank you for your patronage!



Ordering Information

Some products may not be available outside the U.S.A.
Algunos productos no se comercializan fuera de EE.UU.
Certains produits ne sont pas disponibles en dehors des États-Unis.
Einige Produkte sind außerhalb der USA nicht erhältlich.

To Place an Order:

Contact your local BIOMET 3i Representative or call:

BIOMET 3i Customer Service
Monday–Thursday 8am–8pm (EST)
Friday 8am–6:30pm (EST)
800-342-5454, Fax 561-776-1272
In Canada: 800-363-1980, Outside US: 561-776-6700

TERMS AND CONDITIONS OF SALE. The following are the terms and conditions by which BIOMET 3i sells its products.

GOVERNING TERMS. Any shipment of products shall be deemed to be on the terms and conditions stated herein. Any and all terms and conditions submitted by Purchaser are hereby rejected.

We understand that this terms are only applicable if ordered directly in the US and if ordered with local companies, local terms will apply.

TAXES. The prices set forth herein do not include any sales, use, excise, ad valorem, property or other taxes applicable to the sale, use or delivery of the products, all of which shall be paid by Purchaser separately or added to the contract price and paid by Purchaser to BIOMET 3i.

PRICES. The prices set forth for the products are in United States Dollars.

DELIVERY. Products sold hereunder shall be delivered by BIOMET 3i, F.O.B. Palm Beach Gardens, Florida.

PAYMENT. Payment terms shall be as set forth under the Terms section on front of the invoice.

WARRANTY. NEITHER BIOMET 3i NOR ANY OF ITS AFFILIATES MAKE ANY WARRANTIES, EXPRESS OR IMPLIED, EXCEPT THAT ITS PRODUCTS SHALL BE FREE FROM DEFECTS IN MATERIAL AND/OR WORKMANSHIP. This warranty applies only to the original Purchaser. In the event of a product defect, immediately notify BIOMET 3i of the defect prior to returning the product. Devices shall be sterilized prior to return. As its sole obligation and as Purchaser's sole remedy under this warranty, BIOMET 3i will, at its option, either repair, replace or issue credit for such products. Purchaser assumes all risks and liability resulting from the use of these products, whether used separately or in combination with other products.

EXCEPT FOR THE WARRANTY EXPRESSLY DESCRIBED IN THE WARRANTY PARAGRAPH ABOVE, NEITHER BIOMET 3i NOR ANY OF ITS AFFILIATES MAKE ANY OTHER WARRANTY WITH RESPECT TO THE BIOMET 3i PRODUCTS, EXPRESSED OR IMPLIED, WRITTEN OR ORAL, INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR A WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE.

EDUCATION; MODIFICATION OF PRODUCTS. BIOMET 3i reserves the right to modify products, discontinue products and change specifications of products.

RETURN POLICY. A Purchaser may return any BIOMET 3i Product within 90 days of the invoice date. BIOMET 3i Products returned within 60 days of the invoice date will be accepted without any re-stocking fee to the Purchaser. BIOMET 3i Products returned between 61 and 90 days of the invoice date will be subject to a 15% restocking fee. A BIOMET 3i Product will be accepted for return only if the BIOMET 3i Product is (1) returned in its original, unopened package (including autoclavable products); (2) received by BIOMET 3i within 90 days after the invoice date for that product; and (3) returned freight prepaid to 4555 Riverside Drive, Palm Beach Gardens, Florida 33410, (Attention: Returns Department) accompanied by the BIOMET 3i Return Authorization Number provided by the BIOMET 3i Customer

Service Department. Purchaser should be advised that BIOMET 3i Products not meeting the above criteria will not be accepted for return.

EXCHANGE POLICY. BIOMET 3i Products may be exchanged within 180 days of the invoice date, without a re-stocking fee, for other BIOMET 3i Products of equal or lesser value within the same product family as determined by BIOMET 3i (ex. implant for implant, abutment for abutment, membrane for membrane). BIOMET 3i Products may be exchanged within 180 days of the invoice date for other BIOMET 3i Products of greater value within the same product family as determined by BIOMET 3i (ex. implant for implant, abutment for abutment, membrane for membrane), with the price differential paid by the Purchaser. In both cases, the exchanged BIOMET 3i Product must be (1) returned in its original, unopened package (including autoclavable products); (2) a currently offered BIOMET 3i Product (discontinued product not included); and (3) returned freight prepaid to 4555 Riverside Drive, Palm Beach Gardens, Florida 33410, (Attention: Returns Department) and accompanied by the BIOMET 3i Products Return Authorization Number provided by the BIOMET 3i Customer Service Department.

AVAILABILITY OF PRODUCTS. Not all products are registered or available in every region/country. Please contact your BIOMET 3i Representative for product availability and additional information. Offers and products void where prohibited by law.

EXPORT CONTROL. All products are subject to all applicable export control and sanctions legislation.

FORCE MAJEURE. BIOMET 3i shall not be liable for any loss or damage due to failure or delay arising out of any cause beyond the reasonable control of BIOMET 3i. In the event of any failure or delay resulting from such causes, an equitable adjustment of delivery and any other appropriate terms and conditions shall be made. No such failure or delay shall be the basis for an increase of obligation by BIOMET 3i nor any termination by the Purchaser.

LIMITATION OF LIABILITY. Under no circumstances shall BIOMET 3i nor any of its affiliates be liable for business interruption, loss of profits, special, indirect, incidental or consequential damages of any nature and from any cause whether based in contract, warranty or tort (including negligence) or other legal theory, even if BIOMET 3i has been advised of the possibility of such damages.

VALIDITY. If any provision of these Terms and Conditions is found to be illegal or unenforceable in any respect, such illegality or unenforceability shall not affect any other provision of these Terms and Conditions, all of which shall remain enforceable in accordance with their terms.

GENERAL PROVISIONS. The purchase of the products and these Terms and Conditions are governed by the laws of Florida without regard to conflicts of law principles. This document incorporates all oral and written representations between the parties and constitutes the entire agreement and understanding of the parties with respect to the subject matter hereof and supersedes any and all other agreements either oral or written between the parties with respect to such subject matter. No amendment or modification may be made to this document unless made in writing and duly executed by an authorized representative of BIOMET 3i.

CAUTION. Federal (U.S.) Law restricts these devices to sale by or on the order of a licensed healthcare practitioner. Devices labeled STERILE are certified to be sterile unless sterile package is opened or damaged.

EXPORT ORDERS. Require advance payment or letter of credit. Shipments are made either by freight collect or paid by BIOMET 3i and added to the invoice. Unless special instructions are received with an order, the shipping method will be determined by BIOMET 3i, F.O.B. Palm Beach Gardens, Florida.

DUNS #: 186127825
BIOMET 3i is listed in the SAM (U.S. System for Award Management) Database.
©2016 BIOMET 3i, LLC.

Want More? Solutions Start Here!

Check Out The BIOMET **3i** Online Store at
www.shopbiomet3i.com

(Available for France, Spain & US Customers Only)

*For More Information, Please Contact Your Local BIOMET **3i** Sales Representative.*



Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Toll Free: 1-800-342-5454
Phone: +1-561-776-6700
Fax: +1-561-776-1272
www.biomet3i.com



BIOMET **3i** Dental Ibérica S.L.
EMEA Headquarters
WTC Almeda Park, Ed. 4, Planta 2
C/Tirso de Molina, 40
08940, Cornellà de Llobregat (Barcelona)
Spain
Phone: +34-93-470-55-00
Fax: +34-93-371-78-49
www.biomet3i.com



Join
Us



Follow
Us



Watch
Us



Download
It

3i T3, BellaTek, BellaTek design, Certain, DIEM, Encode, EP, GingiHue, Gold-Tite, IOL, Navigator, OSSEOTITE, PreFormance, PREVAIL, Provide, QuickBridge, QuickSeat and ZiReal are registered trademarks of BIOMET **3i**, LLC. NanoTite, Providing Solutions – One Patient At A Time and Twist Lock are trademarks of BIOMET **3i**, LLC. DOLDER is a registered trademark of Dr. Charlotte Fleischer-Dolder and Dr. Fritz Dolder of Switzerland. LOCATOR is a registered trademark of Zest Anchors, Inc. ©2017 BIOMET **3i**, LLC. All rights reserved.

All trademarks herein are the property of BIOMET **3i**, LLC unless otherwise indicated. This material is intended for clinicians only and is NOT intended for patient distribution. This material is not to be redistributed, duplicated or disclosed without the express written consent of BIOMET **3i**. For additional product information, including indications, contraindications, warnings, precautions and potential adverse effects, please visit the BIOMET **3i** Website: www.ifu.biomet3i.com.



INSTRES
REV H 12/17